



## Product Specification Sheet

# VWR MICROSCOPE SLIDES

VWR International certifies that the specified products meet all VWR specifications and quality criteria for release. The products comply in every respect with the relevant particular specifications, drawings, relevant standards and regulations in force.

*VWR International certifie que les produits mentionnés respectent toutes les spécifications et critères de qualité VWR. Les produits satisfont aux spécifications particulières, aux plans, aux normes et règlements en vigueur s'y rapportant.*

VWR International bestätigt hiermit, daß die angegebenen Produkte alle VWR Spezifikationen und Qualitätskriterien zur Freigabe erfüllen. Die Produkte entsprechen in jeder Hinsicht den maßgeblichen Spezifikationen, den Zeichnungen sowie den relevanten Normen und gültigen Vorschriften.

| Cat. No. | Description                | Thickness    | W×D×H (mm) | Pk |
|----------|----------------------------|--------------|------------|----|
| 631-0113 | No frosted end, ISO 8037/1 | 0.8 - 1.0 mm | 76×26×1    | 50 |
| 631-0114 | No frosted end, ISO 8037/1 | 1.0 - 1.2 mm | 76×26×1    | 50 |
| 631-0116 | Frosted end, ISO 8027/1    | 0.8 - 1.0 mm | 76×26×1    | 50 |
| 631-0117 | Frosted end, ISO 8037/1    | 1.0 - 1.2 mm | 76×26×1    | 50 |
| 631-0905 | No frosted end, BS 7011    | 1.0 - 1.2 mm | 76×26×1    | 50 |
| 631-0906 | Frosted end, BS 7011       | 1.0 - 1.2 mm | 76×26×1    | 50 |
| 631-0907 | Frosted end, BS 7011       | 0.8 - 1.0 mm | 76×26×1    | 50 |
| 631-0908 | No frosted end, BS 7011    | 0.8 - 1.0 mm | 76×26×1    | 50 |
| 631-0909 | White frosted end          | 1.0 mm       | 76×26×1    | 50 |

### Storage:

- Keep dry
- Avoid large variations in temperature during both storage and usage. Cooling of the product may cause condensation and lead to condensed water forming between the glasses.
- Slide storage periods: see "use before" date

### Other characteristics:

- Resistant to common solvents
- Twinfrosted slides with embossed surface, for easier writing on and reading
- Frosted surface separates the slides in storage preventing scratching
- These slides are made from low iron content, white glass - ground 90°, washed and polished

**Disclaimer:** VWR wants to state that these specifications can only be regarded as general information and will not discharge the user from his obligation to ensure the product is suitable for the intended use - If there's any doubt ask us.

Date of issue: **November 2012**

VWR International bvba/sprl, Haasrode Research Park Zone 2020, Geldenaaksebaan 464, 3001 Leuven, Belgium

<http://www.vwr.com>



Version: 2.0