

# Certificate of Conformity

|                                         |                                                                                                                             |
|-----------------------------------------|-----------------------------------------------------------------------------------------------------------------------------|
| <b>VWR Catalogue Number</b>             | 631-1160                                                                                                                    |
| <b>Description</b>                      | VWR®, Plain Microscope Slides, Ground edges, clipped corner (45°) slides, 1.0 (-0.05/+0.1) mm thickness, 76×26 mm dimension |
| <b>Lot Number</b>                       | 25/08/2024                                                                                                                  |
| <b>Date of Manufacture (yyyy-mm-dd)</b> | 2024-08-25                                                                                                                  |
| <b>Expiry Date (yyyy-mm-dd)</b>         | 2029-08-25                                                                                                                  |
| <b>Country of Origin</b>                | Manufactured in China                                                                                                       |
| <b>Date of Issue (yyyy-mm-dd)</b>       | 2025-05-21                                                                                                                  |

## Quality System Compliance

Products are manufactured under the **ISO 9001:2015** standard. Products are inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP.

## QC Testing

Representative products samples are collected and inspected in accordance with current applicable QC requirements.

## Product Specifications

|                         |                                                                                                             |
|-------------------------|-------------------------------------------------------------------------------------------------------------|
| <b>Material</b>         | Sheet Glass.                                                                                                |
| <b>Length</b>           | 75.60~ 75.90 mm                                                                                             |
| <b>Width</b>            | 25.00 ~ 25.30 mm                                                                                            |
| <b>Thickness</b>        | 1.00 (-0.05/+0.1) mm                                                                                        |
| <b>Refractive index</b> | 1.53 +/-0.02 of refractive index of the light at the green mercury e-line (wavelength=547.07nm) ISO 8037-2. |

|                                |                                                                                                                                                                                                                                                                                                                                                             |
|--------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Wettability</b>             | The central parts of the two surfaces of the slides should be wetted equably, and slightly dehydration of the edges is allowed. (ISO 8037-2)                                                                                                                                                                                                                |
| <b>Sterilization</b>           | Not applicable.                                                                                                                                                                                                                                                                                                                                             |
| <b>ATP Assay</b>               | Not applicable.                                                                                                                                                                                                                                                                                                                                             |
| <b>DNase &amp; RNase Free</b>  | Not applicable.                                                                                                                                                                                                                                                                                                                                             |
| <b>BSE/TSE</b>                 | No use of any raw material produced from or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE). |
| <b>Cytotoxicity</b>            | Not applicable.                                                                                                                                                                                                                                                                                                                                             |
| <b>Non-Pyrogenic Statement</b> | Not applicable.                                                                                                                                                                                                                                                                                                                                             |
| <b>Latex Statement</b>         | Not applicable.                                                                                                                                                                                                                                                                                                                                             |
| <b>BPA Statement</b>           | Not applicable.                                                                                                                                                                                                                                                                                                                                             |
| <b>DEHP Statement</b>          | Not applicable.                                                                                                                                                                                                                                                                                                                                             |
| <b>RoHS</b>                    | Not applicable.                                                                                                                                                                                                                                                                                                                                             |
| <b>REACH Statement</b>         | Not applicable.                                                                                                                                                                                                                                                                                                                                             |
| <b>Storage Conditions</b>      | Room temperature.                                                                                                                                                                                                                                                                                                                                           |
| <b>Shelf Life</b>              | 5 years.                                                                                                                                                                                                                                                                                                                                                    |

Disclaimer: VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. The purpose of the product is for use in laboratory only.  
This document has been produced electronically and is valid without a signature.