



VWR® 25mL Tissue Culture Flasks

Description: 25mL, Surface Treated, w/plug seal cap or w/vent cap, Sterilized

Purpose: For use in applications where cell attachment is desired


Materials

Flask: GPPS (General polystyrene) Color: Clear

Cap: HDPE(High Density Polyethylene) Color: Green

Membrane: 0.22µm Hydrophobic PP (Only use for Vent Cap)

Dimensions (Unit: mm)



Features

- Temperature: -20°C~+50°C
- Sterilized by irradiation SAL10⁶ (ISO11137)
- DNase/RNase free and Non-pyrogenic
- Bulk packed in zip-Lock bag
- Autoclavability: No
- Maximum working volume: 8 mL
- Recommended working volume: 3 mL
- Growth area: 12.5 cm²
- Every product is printed with a batch number for traceability
- Shelf Life: 3 years after month of production
- Manufactured in a Class 100,000 cleanroom environment
- Manufactured under ISO13485 and ISO9001 quality management system



VWR® 25mL Tissue Culture Flasks

VWR EU Cat. No.	VWR NA Cat. No.	Description	Sterile	Qty/Bag	Qty/Pk	Qty/Case
734-2310	10062-870	25mL Tissue Culture Flask, Surface Treated, w/plug seal cap, Sterilized	Yes	10	50	200
734-2309	10062-868	25mL Tissue Culture Flask, Surface Treated, w/vent cap, Sterilized	Yes	10	50	200



VWR® 50mL Tissue Culture Flasks

Description: 50mL, Surface Treated, w/plug seal cap or w/vent cap, Sterilized

Purpose: For use in applications where cell attachment is desired

Materials

Flask: GPPS (General polystyrene) Color: Clear


Cap: HDPE(High Density Polyethylene) Color: Green

Membrane: 0.22µm Hydrophobic PP (Only use for Vent Cap)

Features

- Temperature: -20°C~+50°C
- Sterilized by irradiation SAL10⁶ (ISO11137)
- DNase/RNase free and Non-pyrogenic
- Bulk packed in zip-Lock bag
- Autoclavability: No
- Maximum working volume: 17.5 mL
- Recommended working volume: 6.25 mL
- Growth area: 25.0 cm²
- Every product is printed with a batch number for traceability
- Shelf Life: 3 years after month of production
- Manufactured in a Class 100,000 cleanroom environment
- Manufactured under ISO13485 and ISO9001 quality management system

Dimensions (Unit: mm)



STERILE

R

VWR® 50mL Tissue Culture Flasks

VWR EU Cat. No.	VWR NA Cat. No.	Description	Sterile	Qty/Bag	Qty/Pk	Qty/Case
734-2312	10062-874	50mL Tissue Culture Flask, Surface Treated, w/plug seal cap, Sterilized	Yes	10	50	200
734-2311	10062-872	50mL Tissue Culture Flask, Surface Treated, w/vent cap, Sterilized	Yes	10	50	200



VWR® 250mL Tissue Culture Flasks

Description: 250mL, Surface Treated, w/plug seal cap or w/vent cap, Sterilized

Purpose: For use in applications where cell attachment is desired

Materials

Flask: GPPS (General polystyrene) Color: Clear


Cap: HDPE(High Density Polyethylene) Color: Green

Membrane: 0.22µm Hydrophobic PP (Only use for Vent Cap)

Features

- Temperature: -20°C~+50°C
- Sterilized by irradiation SAL10⁻⁶ (ISO11137)
- DNase/RNase free and Non-pyrogenic
- Bulk packed in zip-Lock bag
- Autoclavability: No
- Maximum working volume: 60 mL
- Recommended working volume: 18.75 mL
- Growth area: 75.0 cm²
- Every product is printed with a batch number for traceability
- Shelf Life: 3 years after month of production
- Manufactured in a Class 100,000 cleanroom environment
- Manufactured under ISO13485 and ISO9001 quality management system

Dimensions (Unit: mm)



STERILE

R

VWR® 250mL Tissue Culture Flasks

VWR EU Cat. No.	VWR NA Cat. No.	Description	Sterile	Qty/Bag	Qty/Pk	Qty/Case
734-2314	10062-862	250mL Tissue Culture Flask, Surface Treated, w/plug seal cap, Sterilized	Yes	5	25	100
734-2313	10062-860	250mL Tissue Culture Flask, Surface Treated, w/vent cap, Sterilized	Yes	5	25	100



VWR® 600mL Tissue Culture Flasks

Description: 600mL, Surface Treated, w/plug seal cap or w/vent cap, Sterilized

Purpose: For use in applications where cell attachment is desired

Materials

Flask: GPPS (General polystyrene) Color: Clear


Cap: HDPE(High Density Polyethylene) Color: Green

Membrane: 0.22µm Hydrophobic PP (Only use for Vent Cap)

Features

- Temperature: -20°C~+50°C
- Sterilized by irradiation SAL10⁶ (ISO11137)
- DNase/RNase free and Non-pyrogenic
- Bulk packed in zip-Lock bag
- Autoclavability: No
- Maximum working volume: 125 mL
- Recommended working volume: 45.5 mL
- Growth area: 182.0 cm²
- Every product is printed with a batch number for traceability
- Shelf Life: 3 years after month of production
- Manufactured in a Class 100,000 cleanroom environment
- Manufactured under ISO13485 and ISO9001 quality management system

Dimensions (Unit: mm)



STERILE

R

VWR® 600mL Tissue Culture Flasks

VWR EU Cat. No.	VWR NA Cat. No.	Description	Sterile	Qty/Bag	Qty/Pk	Qty/Case
734-2316	10062-866	600mL Tissue Culture Flask, Surface Treated, w/plug seal cap, Sterilized	Yes	5	20	40
734-2315	10062-864	600mL Tissue Culture Flask, Surface Treated, w/vent cap, Sterilized	Yes	5	20	40



VWR® 600mL Extended Tissue Culture Flasks

Description: 600mL, extended, Surface Treated, w/plug seal cap or w/vent cap, Sterilized

Purpose: For use in applications where cell attachment is desired


Materials

Flask: GPPS (General polystyrene) Color: Clear

Cap: HDPE(High Density Polyethylene) Color: Green

Membrane: 0.22µm Hydrophobic PP (Only use for Vent Cap)

Dimensions (Unit: mm)



Features

- Extended height reduces the risk of overflowing media
- Temperature: -20°C~+50°C
- Sterilized by irradiation SAL10⁻⁶ (ISO11137)
- DNase/RNase free and Non-pyrogenic
- Bulk packed in zip-Lock bag
- Autoclavability: No
- Maximum working volume: 200 mL
- Recommended working volume: 45.5mL
- Growth area: 182.0 cm²
- Meet the requirements of BSE/TSE
- Latex Statement: this product is latex free
- Every product is printed with a batch number for traceability
- Shelf Life: 3 years after month of production
- Manufactured in a Class 100,000 cleanroom environment
- Manufactured under ISO13485 and ISO9001 quality management system



STERILE

R

VWR® 600mL Extended Tissue Culture Flasks

VWR EU Cat. No.	VWR NA Cat. No.	Description	Sterile	Qty/PK	Qty/Box	Qty/Case
734-3265		600mL Tissue Culture Flask, extended, Surface Treated, w/plug seal cap, Sterilized	Yes	5	20	40
734-3266		600mL Tissue Culture Flask, extended, Surface Treated, w/vent cap, Sterilized	Yes	5	20	40



VWR® 850mL Tissue Culture Flasks

Description: 850mL, Surface Treated, w/plug seal cap or w/vent cap, Sterilized

Purpose: For use in applications where cell attachment is desired

Materials

Flask: GPPS (General polystyrene) Color: Clear


Cap: HDPE(High Density Polyethylene) Color: Green

Membrane: 0.22µm Hydrophobic PP (Only use for Vent Cap)

Features

- Temperature: -20°C~+50°C
- Sterilized by irradiation SAL10⁶ (ISO11137)
- DNase/RNase free and Non-pyrogenic
- Bulk packed in zip-Lock bag
- Autoclavability: No
- Maximum working volume: 200 mL
- Recommended working volume: 75 mL
- Growth area: 300.0 cm²
- Every product is printed with a batch number for traceability
- Shelf Life: 3 years after month of production
- Manufactured in a Class 100,000 cleanroom environment
- Manufactured under ISO13485 and ISO9001 quality management system

Dimensions (Unit: mm)



VWR® 850mL Tissue Culture Flasks

VWR EU Cat. No.	VWR NA Cat. No.	Description	Sterile	Qty/Bag	Qty/Pk	Qty/Case
734-2601	10062-886	850mL Tissue Culture Flask, Surface Treated, w/plug seal cap, Sterilized	Yes	3	9	18
734-2600	10062-884	850mL Tissue Culture Flask, Surface Treated, w/vent cap, Sterilized	Yes	3	9	18