

# Certificate of Conformity

|                                   |                                                                                                                                                                                               |
|-----------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>VWR Catalogue Number</b>       | 662-1792                                                                                                                                                                                      |
| <b>Description</b>                | Speciality pH electrodes, with 1 m fixed cable for Water, creams, beverage, pH range: 0 to 14, pH only half cell type, 0...+100 °C temperature, 12 mm diameter, 150 mm length, DIN connection |
| <b>Lot Number</b>                 | 4523                                                                                                                                                                                          |
| <b>Country of Origin</b>          | Manufactured in USA                                                                                                                                                                           |
| <b>Date of Issue (yyyy-mm-dd)</b> | 2024-01-09                                                                                                                                                                                    |

## Quality System Compliance

Products are manufactured under the **ISO 9001:2015** standards. Products are inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP.

## QC Testing

Representative product sensors are individually tested according to quality procedure DOP-Q126.

## Product Specifications

|                                |                                                                                                                |
|--------------------------------|----------------------------------------------------------------------------------------------------------------|
| <b>Material</b>                | Glass.                                                                                                         |
| <b>Sterilization</b>           | Product is chemical treated and they meet a minimum requirement of 10 <sup>-6</sup> Sterility assurance level. |
| <b>ATP Assay</b>               | Not Applicable.                                                                                                |
| <b>DNase &amp; RNase Free</b>  | Not Applicable.                                                                                                |
| <b>BSE/TSE</b>                 | Not Applicable.                                                                                                |
| <b>Cytotoxicity</b>            | Not Applicable.                                                                                                |
| <b>Non-Pyrogenic Statement</b> | Not Applicable.                                                                                                |
| <b>Latex Statement</b>         | Not Applicable.                                                                                                |

|                           |                                                                                                                                                                                                                                                                                                                                                                      |
|---------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>BPA Statement</b>      | Not Applicable.                                                                                                                                                                                                                                                                                                                                                      |
| <b>DEHP Statement</b>     | Please see RoHS information below.                                                                                                                                                                                                                                                                                                                                   |
| <b>RoHS</b>               | No substances (Lead, Cadmium, Mercury, Hexavalent Chromium (Cr6+), Poly Brominated Biphenyls (PBB), Poly Brominated Diphenyl ethers (PBDE), Bis(2-Ethyhexyl) Phthalate (DEHP), Benzyl Butyl Phthalate (BBP), Dibutyl Phthalate (DBP), Diisobutyl Phthalate (DIBP)) are used in manufacturing the raw materials and final product. No routinely analyse is performed. |
| <b>REACH Statement</b>    | No substances as listed on >Substance of Very High Concern< candidate list published 2023-06-14 are used in manufacturing of the raw materials used to produce this product. No routinely analysis is performed.                                                                                                                                                     |
| <b>Storage Conditions</b> | 0° C - 50° C, wet.                                                                                                                                                                                                                                                                                                                                                   |

Disclaimer: VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. The purpose of the product is for use in laboratory only.  
This document has been produced electronically and is valid without a signature.