

Certificate of Conformity

VWR Catalogue Number	See List
Description	Plastic Box, Translucent Polypropylene
Country of Origin	Manufactured in France
Date of Issue (yyyy-mm-dd)	2025-11-20

VWR Catalogue Number	W x D x H ext [mm]	Description	Pk
216-3024	90 × 60 × 20	Plastic box, 1 compartment, compartment size 82×53 mm	24
216-3025	90 × 60 × 20	Plastic box, 2 compartments, compartment size 53×40 mm	24
216-3026	90 × 60 × 20	Plastic box, 4 compartments, compartment size 40×26 mm	24
216-3027	90 × 60 × 20	Plastic box, 6 compartments, compartment size 26×26 mm	24
216-3028	130 × 90 × 25	Plastic box, 1 compartment, compartment size 121×72 mm	16
216-3029	130 × 90 × 25	Plastic box, 12 compartments, compartment size 29×26 mm	16
216-3030	190 × 130 × 30	Plastic box, 1 compartment, compartment size 171×118 mm	16
216-3031	190 × 130 × 30	Plastic box, 12 compartments, compartment size 40×38 mm	16
216-3032	260 × 180 × 40	Plastic box, 1 compartment, compartment size 250×153 mm	14
216-3033	260 × 180 × 40	Plastic box, 24 compartments, compartment size 39×38 mm	14
216-3034	360 × 260 × 50	Plastic box, 1 compartment, compartment size 351×239 mm	6
216-3035	360 × 260 × 50	Plastic box, 28 compartments, compartment size 59×45 mm	6

Quality System Compliance

This certificate confirms that the product has been sourced from a supplier or manufacturer who has declared their Quality Management System (QMS) complies with **ISO 14001:2015**

QC Testing

Representative products samples are collected and inspected in accordance with current applicable QC requirements.

Product Specifications

Material Polypropylene.

Sterilization	Non-Sterile.
ATP Assay	Not Applicable.
DNase & RNase Free	Not Applicable.
BSE/TSE	No use of any raw material produced from or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE).
Cytotoxicity	Not Applicable.
Non-Pyrogenic Statement	Not Applicable.
Latex Statement	Does not intentionally incorporate in the composition.
BPA Statement	Does not intentionally incorporate in the composition of its polymers.
DEHP Statement	This product does not contain intentionally in its composition any phthalate.
RoHS	Raw materials comply with the requirements of the Directive 2011/65/EU which establishes a maximum concentration value upto 0.1% by weight for lead, mercury, hexavalent chromium (VI), Polybrominated Biphenyls (PBB) and Polybrominated Diphenyl Ethers (PBDE), and up to 0.01% for cadmium.
REACH Statement	Product meets the criteria for authorization under the European Union's REACH regulation (EC) No. 1907/2006. To the best of our knowledge product does not contain SVHC's above the regulatory threshold of 0.1%.
Storage Conditions	< 60°.
Shelf Life	10 years.

Disclaimer: VWR's Quality Management System (QMS) is ISO 9001:2015 certified as a wholesale distributor of finished products, ensuring effective oversight of finished product handling, traceability, and distribution processes. VWR is not a manufacturer and does not engage in manufacturing, co-manufacturing or processing activities such as design, packaging, labeling, testing or quality release. Consequently, the supplier or manufacturer's QMS certification may differ from VWR's.

VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. This product is intended for research use only. It is not for use in diagnostic procedures, therapeutic applications, or any human or clinical use.

This document has been produced electronically and is valid without a signature.

Version History

Sl.no	Version number	Description	Month & Year
1	1.0	Created General CoC using new layout	11-2025