

Sodium Hydroxide
Baker
Pellets



Material No: 0404.1000
Effective from: 2013/04/17

Product Specification

Test	Specification
Assay	$\geq 97 \%$
Aluminium (Al)	$\leq 0.002 \%$
Carbonate (as Na_2CO_3)	$\leq 1 \%$
Chloride (Cl)	$\leq 0.01 \%$
Heavy Metals (as Pb)	$\leq 0.002 \%$
Iron (Fe)	$\leq 0.002 \%$
Sulphate (SO_4)	$\leq 0.01 \%$

For Laboratory, Research or Manufacturing use.

Barbara Rykowska
Quality Assurance Coordinator

For questions on this Product Specification please contact Technical Service at technical.service@avantorsciences.com
Avantor Performance Materials Poland S.A. Sowinskiego 11 Gliwice 44-101 Poland Tel: +48 322392000