

Certificate of Conformity

VWR Catalogue Number	See List
Description	Stirring shafts, VWR®
Country of Origin	Manufactured in United Kingdom
Date of Issue (yyyy-mm-dd)	2022-08-12

Catalogue Number	Length [mm]	Diameter [mm]
441-0253	300	8
441-0254	400	8
441-0255	500	8
441-0256	650	10
441-0257	750	10
441-0258	1000	10
441-0259	650	12
441-0260	750	12
441-0261	1000	12

Quality System Compliance

Products are manufactured under the **ISO 9001:2015** standard. Products are inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable QC requirements.

Product Specifications

Material

Virgin Polytetrafluoroethylene (PTFE) Grade TF1620.
Raw Material complies with Directive FDA 21 CFR 177.1550 and complies with USP Class VI - USP 27, NF 22,2004 for plastic at 70 °C.
The test is carried out by 3M-Dyneon who manufactures PTFE used and refers to PTFE raw material.

Temperature Resistance	-200 °C to +280 °C.
Sterilization	Non-Sterile.
ATP Assay	Not Applicable.
DNase & RNase Free	Not Applicable.
BSE/TSE	With reference to guideline EMEA/410/01.rev.3 the above product contains no ingredients of animal or human origin and does not come in contact with material derived from or exposed to animals affected by or under quarantine for Transmitting, Animal Spongiform Encephalopathy (TSE)/ Bovine Spongiform Encephalopathy (BSE).
Cytotoxicity	Not Applicable.
Non-Pyrogenic Statement	Not Applicable.
Latex Statement	This product is not made from natural rubber latex.
BPA Statement	Bisphenol are not used in the manufacture of the raw material and are not expected to be present.
DEHP Statement	Please see RoHS information below.
RoHS	No substances (Lead, Cadmium, Mercury, Hexavalent Chromium (Cr6+), Poly Brominated Biphenyls (PBB), Poly Brominated Diphenyl ethers (PBDE), Bis(2-Ethyhexyl) Phthalate (DEHP), Benzyl Butyl Phthalate (BBP), Dibutyl Phthalate (DBP), Diisobutyl Phthalate (DIBP)) are used in manufacturing the raw materials and final product. No routinely analyse is performed.
REACH Statement	No substances as listed on >Substance of Very High Concern< candidate list published 2022-01-17 are used in manufacturing of the raw materials used to produce this product. No routinely analysis is performed.
Storage Conditions	Ambient Temperature.
Shelf Life	Indefinite when stored under ideal conditions.

Disclaimer: VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. The purpose of the product is for use in laboratory only.
 This document has been produced electronically and is valid without a signature.