

Certificate of Conformity

VWR Catalogue Number	See List
Description	Adhesive PCR Foil Sheet, 63.5 µm Thickness
Country of Origin	Manufactured in USA
Date of Issue (yyyy-mm-dd)	2021-12-02

Catalogue Number	Sterile	Optical Property
391-1257	Sterile	Light blocking
391-1256	Non-Sterile	Light blocking

Quality System Compliance

Products are inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable QC requirements.

Product Specifications

Material Aluminium Foil.

Recommended Temperature -80 °C to 120 °C.

Sterilization Product labelled as sterile is gamma irradiated and meets SAL 10⁻³ with a required dose range of 12-25 kGy. (which is not applicable for Non-sterile products)

ATP Assay Not applicable.

DNase & RNase Free This product is free of any detectable DNase/RNase contamination, and Protease.

DNA Contamination This product is free of detectable levels of Human/Mouse DNA contamination (PCR Certified).

BSE/TSE	No use of any raw material produced from or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE).
Cytotoxicity	Testing is conducted and have been shown to be non-toxic.
Non-Pyrogenic Statement	Product is non-pyrogenic.
Latex Statement	This product is not made from natural rubber latex.
BPA Statement	Not applicable.
DEHP Statement	Not applicable.
RoHS	No substances (Lead, Cadmium, Mercury, Hexavalent Chromium (Cr ⁶⁺), Poly Brominated Biphenyls (PBB), Poly Brominated Diphenyl ethers (PBDE), Bis(2-Ethyhexyl) Phthalate (DEHP), Benzyl Butyl Phthalate (BBP), Dibutyl Phthalate (DBP), Diisobutyl Phthalate (DIBP)) are used in manufacturing the raw materials and final product. No routinely analyse is performed.
REACH Statement	No substances as listed on >Substance of Very High Concern< candidate list published 2021-07-08 are used in manufacturing of the raw materials used to produce this product. No routinely analyse is performed.
Storage Conditions	21.1 °C up to -12 °C; 50% RH ± 10%
Shelf Life	3 years

Disclaimer: VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. The purpose of the product is for use in laboratory only.
This document has been produced electronically and is valid without a signature.