

Certificate of Conformity

VWR Catalogue Number	See List
Description	VWR®, Reagent Bottles with PP 28 Thread
Country of Origin	Manufactured in Germany
Date of Issue (yyyy-mm-dd)	2024-01-18

Clear Glass

Catalogue Number	Capacity [mL]	Thread	Diameter Ext [mm]	Height [mm]
215-4403	30	PP 28	34.2	69
215-4404	50	PP 28	40.5	81
215-4405	100	PP 28	50.5	96.5
215-4406	200	PP 28	62.5	119.5
215-4407	250	PP 28	65	133
215-4408	500	PP 28	82	161
215-4409	1000	PP 28	100.5	206.2

Amber Glass

Catalogue Number	Capacity [mL]	Thread	Diameter Ext [mm]	Height [mm]
VWRI215-4410_P	30	PP 28	34.2	69
215-4411	50	PP 28	40.5	81
215-4412	60	PP 28	38.6	82.7
215-4413	100	PP 28	50.5	96
215-4414	150	PP 28	56.7	111
215-4415	200	PP 28	62.5	119.5
215-4416	250	PP 28	65	133
215-4417	500	PP 28	82	161
215-4418	1000	PP 28	100.5	106.2

Screw caps for PP 28 bottles

Catalogue Number	Description	Thread	Colour
215-4458	Screw cap, PP, with foam insert	PP 28	Black
215-4459	Screw cap, PP, with foam insert	PP 28	White

Quality System Compliance

Products are manufactured under the **DIN EN ISO 9001:2015** standard. Products are inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable QC requirements.

Product Specifications

Material	Bottle: Clear or Amber soda-lime glass. Screw Cap: Polypropylene.
Sterilization	Non-Sterile.
FDA	The materials used to produce the products are listed in List of Authorised Substances of Regulation EU/10/2011. The material grades conform to EU Directive 1935/2004 and 2023/2006 on materials and articles intended to come into contact with food.
ATP Assay	Not applicable.
DNase & RNase Free	Not applicable.
BSE/TSE	No use of any raw material produced from or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE).
Cytotoxicity	Not Applicable.
Non-Pyrogenic Statement	Not applicable.

Latex Statement	The above products are not made from natural rubber latex.
BPA Statement	Bisphenol are not used in the manufacture of the raw material and are not expected to be present.
DEHP Statement	Not applicable.
RoHS	Not applicable.
REACH Statement	No substances as listed on >Substance of Very High Concern< candidate list published 2023-01-17 are used in manufacturing of the raw materials used to produce this product. No routinely analysis is performed.
Storage Conditions	No special conditions.
Shelf Life	Not ascertainable, too long.

Disclaimer: VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. The purpose of the product is for use in laboratory only.
This document has been produced electronically and is valid without a signature.

Version History

Sl.no	Version number	Description	Month & Year
1	1.0	Created General CoQ using new layout	07-2023
2	1.1	Added FDA information	01-2024