

Certificate of Quality

VWR Catalogue Number	See List
Description	VWR® Culture Tubes, Plastic, without Caps
Country of Origin	Manufactured in USA
Date of Issue (yyyy-mm-dd)	2026-06-09

Polypropylene Material

VWR Catalogue Number	Capacity	Height	Diameter	Packaging	Sterility	Unit
60818-430	5 mL	75 mm	12 mm	Bulk	Non sterile	Case of 1,000
60818-281	5 mL	75 mm	12 mm	125/Bag	Non sterile	Case of 1,000
60818-306	5 mL	75 mm	12 mm	125/Bag	Sterile	Case of 1,000
60818-860	8 mL	100 mm	13 mm	Bulk	Non sterile	Case of 1,000
60818-918	8 mL	100 mm	13 mm	125/Bag	Non sterile	Case of 500
60818-930	8 mL	100 mm	13 mm	125/Bag	Sterile	Case of 1,000
60818-618	14 mL	100 mm	17 mm	Bulk	Non sterile	Case of 1,000
60818-626	14 mL	100 mm	17 mm	125/Bag	Non sterile	Case of 1,000
60818-645	14 mL	100 mm	17 mm	125/Bag	Sterile	Case of 1,000

Polystyrene Material

VWR Catalogue Number	Capacity	Height	Diameter	Packaging	Sterility	Unit
60818-408	5 mL	75 mm	12 mm	Bulk	Non sterile	Case of 1,000
60818-361	5 mL	75 mm	12 mm	Bulk	Non sterile	Case of 3,000
60818-270	5 mL	75 mm	12 mm	125/Bag	Non Sterile	Case of 1,000
60818-292	5 mL	75 mm	12 mm	125/Bag	Sterile	Case of 1,000
60818-849	8 mL	100 mm	13 mm	Bulk	Non sterile	Case of 1,000
60818-907	8 mL	100 mm	13 mm	125/Bag	Non sterile	Case of 1,000
60818-929	8 mL	100 mm	13 mm	125/Bag	Sterile	Case of 1,000

60818-615	14 mL	100 mm	17 mm	Bulk	Non sterile	Case of 1,000
60818-614	14 mL	100 mm	17 mm	125/Bag	Non sterile	Case of 1,000
60818-623	14 mL	100 mm	17 mm	125/Bag	Sterile	Case of 1,000

Quality System Compliance

This certificate confirms that the product has been sourced from a supplier or manufacturer who has declared their Quality Management System (QMS) complies with **ISO 9001:2015** standard.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable product specifications.

Product Specifications

Material	Refer table above. Raw materials Complies with FDA 21 CFR and USP class <88> and <87>.
Sterilization	Products labelled as sterile are beta irradiated and meets a minimum requirement of 10 ⁻⁶ SAL with a specified dose range of 15.2 – 32.0 kGy.
ATP Assay	The products are free of ATP.
DNase & RNase Free	The products are free of any detectable DNase/RNase contamination
BSE/TSE	No use of any raw material produced from or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE).
Cytotoxicity	The products are free of cytotoxicity effects.
Non-Pyrogenic Statement	The products are free of pyrogens.
Latex Statement	The products are not made from natural rubber latex.
BPA Statement	Bisphenols are not used in the manufacture of the raw material and are not expected to be present.
DEHP Statement	DEHP is not used in the manufacture of the raw material and are not expected to be present.

RoHS	No substances (Lead, Cadmium, Mercury, Hexavalent Chromium (Cr ⁶⁺), Poly Brominated Biphenyls (PBB), Poly Brominated Diphenyl ethers (PBDE), Benzyl Butyl Phthalate (BBP), Dibutyl Phthalate (DBP), Diisobutyl Phthalate (DIBP)) are used in manufacturing the raw materials and final product. No routine analysis is performed.
REACH Statement	No substances as listed on >Substance of Very High Concern< candidate list are used in manufacturing of the raw materials used to produce this product. No routine analysis is performed.
Storage Conditions	Storage at ambient temperature is recommended.
Shelf-Life	60 months for sterile products. 72 months for non-sterile products.

Disclaimer: VWR'S Quality Management System (QMS) is ISO 9001:2015 certified as a wholesale distributor of finished products, ensuring effective oversight of finished product handling, traceability and distribution processes. VWR is not a manufacturer and does not engage in manufacturing, co- manufacturing, or processing activities such as design, packaging, labelling, testing, or quality release. Consequently, the supplier or manufacturer's QMS certification may differ from VWR's.

VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. This product is intended for research use only. It is not for use in diagnostic procedures, therapeutic applications, or any human or clinical use. This document has been produced electronically and is valid without a signature.

Version History

Sl.no	Version number	Description	Month & Year
1	1.0	Created General CoQ using new layout	June 2026