

3-part Syringe

Product Information Card



Description:

Sterile 3-part syringe for single use of *Terumo (Philippines) Corporation origin*.

The syringe is intended for the injection of fluids - immediately after filling - into the human body, for aspiration of fluids or for the sampling of blood or other body fluids.

Does not contain components made of PVC or natural rubber latex.

Components:

Description	Raw Materials
Barrel	Polypropylene
Plunger	Polystyrene (for 1ml syringe only) Polypropylene
Gasket	Thermoplastic elastomer
Lubricant	Silicone oil



Product Range:

Product Code	Size	Tip	Replacement for
SS+01T1	1 ml	Concentric, Luer Slip	BS-01T
SS+03L1	3 ml	Concentric, Luer Lock	
SS+10ES1	10 ml	Eccentric, Luer Slip	SS-10ES
SS+20ES1	20 ml	Eccentric, Luer Slip	
SS+20L1	20 ml	Concentric, Luer Lock	
SS+50ES1	50 ml	Eccentric, Luer Slip	BS-50ES
SS+50L1	50 ml	Concentric, Luer Lock	BS-50LG
SS+50C1	50 ml	Catheter Tip	BS-50C

Technical Data:

Product Code	Nominal Capacity	Graduated Capacity	Increments
SS+01T1	1 ml	1 ml	0.01
SS+03L1	3 ml	3 ml	0.1
SS+10ES1	10 ml	10 ml	0.2
SS+20ES1	20 ml	20 ml	1.0
SS+20L1	20 ml	20 ml	1.0
SS+50ES1	50 ml	60 ml	1.0
SS+50L1	50 ml	60 ml	1.0
SS+50C1	50 ml	60 ml	1.0

3-Part Syringe

Packaging:

Single unit: blister pack (film/film)

Product Code	Unit Box	Shipping Carton		
		Quantity	Gross Weight	Dimensions (LxWxH)
SS+01T1	100	1800	7.50 kg	430 x 418 x 306 mm
SS+03L1	100	1800	7.50 kg	640 x 385 x 262 mm
SS+10ES1	100	1200	11.00 kg	546 x 454 x 428 mm
SS+20ES1	50	600	10.15 kg	570 x 400 x 366 mm
SS+20L1	50	600	10.15 kg	570 x 400 x 366 mm
SS+50ES1	25	100	5.10 kg	466 x 312 x 268 mm
SS+50L1	25	100	5.10 kg	466 x 312 x 268 mm
SS+50C1	25	100	5.40 kg	466 x 312 x 268 mm

Sterilization:

With e-beam radiation to a Sterility Assurance Level of at least 10^{-6}

Storage:

Shelf life of 5 years

* The batch number and expiry date are embossed in the film of the blister pack

* Do not store at extreme temperature and humidity

Precautions for Use:

Check integrity of the unit package before use

Do not reuse

Dispose of safely after single use to avoid risk of infection

Status of Regulation:

Class IIa, sterile product following MDD 93/42/EEC

Notified body:

TÜV Rheinland 0197

Manufacturer:

TERUMO (PHILIPPINES) CORPORATION

124 East Main Avenue

Laguna Technopark, Binan

Philippines

EC representative:

TERUMO EUROPE NV

Interleuvenlaan 40

B - 3001 Leuven

Belgium

HO 180GB-10.07TE-1(07.16)E



Terumo Corporation
+81 3 3374 8111

Terumo Europe NV
+32 16 38 12 11

Terumo Medical Products EMEA
+32 16 38 12 11

EMEA Sales Offices

**Terumo Europe NV
African Business Division**
+32 16 38 13 08

**Terumo Europe NV
Benelux Sales Division
Belgium**

+32 16 339 25 80
The Netherlands
0800 0220396

**Terumo Europe NV
Emerging Markets Division**
+32 16 38 12 11

Terumo Deutschland GmbH
+49 6196 80 230

Terumo Deutschland GmbH
Switzerland
+41 56 419 10 10

Terumo Europe España SL
+34 902 10 12 98

Terumo France S.A.S.
+33 130 96 13 00

Terumo Italia S.R.L.
+39 0651 96 14 20

Terumo Russia LLC
+7 495 988 4740

Terumo Sweden AB
+46 3174 85 880

**Terumo Sweden AB
Denmark**
+45 7020 93 80

Terumo Middle East FZE
+971 4 292 0200

Terumo UK Ltd
+44 1276 480 440

**Terumo BCT Tibbi Cihazlar Dagitim
ve Hizmetleri A.S.**
+90 216 590 0074

