

Material No.: 75838-228

Material Description: VWR Pro style step pipettor tips, sterile, 50mL, 25/pack

Lot No.: 1019299

Compliance and Procedure Result

Quality System Compliance: Products are manufactured in **ISO 13485:2016** and **ISO 9001:2015** certified plant. Products are inspected and controlled in accordance with the current applicable product specification.

QC Testing: Representative products samples are collected and inspected in accordance with current applicable QC requirements.

Material:

Tube: Polypropylene

Plunger 0.1 ml/ 0.2ml: Polybutylene terephthalate (PBT)

Plunger 0.5 ml - 50.0 ml: high density polyethylene (HDPE).

Sterilization: Sterilization by irradiation under controlled conditions.

ATP Assay: Tested with detection level < 1×10^{-12} mmol

PCR Inhibitor: Product yield equivalent to control < 10 amplifiable targets

DNase & RNase free: Tested with detection level < 1×10^{-6} KunitzU (DNase) and < 1×10^{-9} KunitzU (RNase)

Nucleic Acid Free: Not applicable.

Manufacturing Environment: Dispenser tips are manufactured under controlled conditions.

BSE/TSE: Tallow derived additives may be used in the manufacture of the raw materials. These tallow derived additives are in compliance with the guidance EMA/410/01 Rev. 3 - July 1, 2011 "Minimizing the risk of transmitting animal spongiform encephalopathy agents via human and veterinary medicinal products", adopted by the European Commission and published in the Official Journal of the European Union March 5, 2011 (2011/C 73/01).

To the best of our knowledge, these treatments ensure inactivation of any agents present in the tallow derivatives to be used in plastics materials. We therefore can state that this grade is considered safe to use with respect to BSE and TSE transmissions.

Cytotoxicity: Not applicable

Endotoxin Free: Tested with detection level < 0,02 EU/Unit

Latex Statement: Not applicable

BPA Statement : Not applicable

DEHP Statement: Not applicable

RoHS: The raw materials used are compliant with the currently valid directive 2011/65/EU and 2002/95/EG (RoHS).

REACH Statement: The product does not contain any "SVHC" ingredients, in accordance with the currently valid regulation (EC) No 1907/2006 on REACH

Storage Conditions: Keep dry and away from sunlight.



Certificate of Compliance

5/19/2023

Signed on behalf of VWR By: James Osborn

Quality Employee

5/19/2023