

Certificate of Quality

VWR Catalogue Number	See List
Description	VWR® Weighing Paper
Country of Origin	Manufactured in USA
Date of Issue (yyyy-mm-dd)	2026-03-20

VWR Catalogue Number	Size	Unit
12578-121	7.6 × 7.6 cm (3 × 3")	Pack of 500
		Case of 24000
12578-165	10.2 × 10.2 cm (4 × 4")	Pack of 500
		Case of 36000
12578-201	15.2 × 15.2 cm (6 × 6")	Pack of 500
		Case of 18000

Quality System Compliance

This certificate confirms that the product has been sourced from a supplier or manufacturer who has declared the products are manufactured in accordance with their internal quality system. Products are inspected and controlled through whole production processing in accordance with their applicable product specifications and QC SOP.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable product specifications i.e. Basis Weight, Thickness, Moisture Content, Tear MD / Tear CD, Tensile MD, Opacity, Gloss, pH Level, Porosity.

Product Specifications

Material	25# pHneutral Glassine Paper The above products comply with 21 CFR 176, 170 – Table 2, conditions A-H, 180,175, 178, 181, 186.1.
Sterilization	The products are non-sterile.
ATP Assay	Not Applicable.

DNase & RNase Free	Not Applicable.
BSE/TSE	Not Applicable.
Cytotoxicity	Not Applicable.
Non-Pyrogenic Statement	Not Applicable.
Latex Statement	The products are free from latex.
BPA Statement	Bisphenols are not used in the manufacture of the raw material and are not expected to be present.
DEHP Statement	Not Applicable.
RoHS	No substances (Lead, Cadmium, Mercury, Hexavalent Chromium (Cr6+), Poly Brominated Biphenyls (PBB), Poly Brominated Diphenyl ethers (PBDE), Benzyl Butyl Phthalate (BBP), Dibutyl Phthalate (DBP), Diisobutyl Phthalate (DIBP) are used in manufacturing the raw materials and final product. No routine analysis is performed.
REACH Statement	No substances as listed on >Substance of Very High Concern< candidate list are used in manufacturing of the raw materials used to produce this product. No routine analysis is performed.
Storage Conditions	Store in a cool, dry area away from moisture.
Shelf Life	Indefinite shelf life (if properly stored).

Disclaimer: VWR'S Quality Management System (QMS) is ISO 9001:2015 certified as a wholesale distributor of finished products, ensuring effective oversight of finished product handling, traceability and distribution processes. VWR is not a manufacturer and does not engage in manufacturing, co- manufacturing, or processing activities such as design, packaging, labelling, testing, or quality release. Consequently, the supplier or manufacturer's QMS certification may differ from VWR's.

VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. This product is intended for research use only. It is not for use in diagnostic procedures, therapeutic applications, or any human or clinical use. This document has been produced electronically and is valid without a signature.

Version History

Sl.no	Version number	Description	Month & Year
1	1.0	Created General CoQ using new layout	March 2026