

Certificate of Quality

VWR Catalogue Number	See List
Description	VWR® Disposable Transfer Pipettes with Reference Lines
Country of Origin	Manufactured in China
Date of Issue (yyyy-mm-dd)	2026-04-21

Capacity 7 ml, Overall length 15.2 cm (6"), Graduations 1 ml, Max. graduation volume 2 ml

VWR Catalog number	Sterility	Packed	Unit
414004-036	Non sterile	500/box, 1000/case	Case of 1,000
			Pallet of 48000
414004-037	Sterile	1/bag, 500/box	Case of 500
			Pallet of 24000

Quality System Compliance

This certificate confirms that the product has been sourced from a supplier or manufacturer who has declared their Quality Management System (QMS) complies with the [ISO13485:2016](#) standard.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable product specifications.

Product Specifications

Material	Low-Density Polyethylene (LDPE) DJ 200A (Virgin)
Sterilization	Product labelled as sterile is EB (Electron beam) irradiated and dose released upon ISO11137 recommended practices in effect at the time of validation. Product labelled as sterile meet a minimum requirement of 10 ⁻⁶ SAL with a specified dose range of 15 – 40 kGy.
ATP Assay	Not Applicable.
DNase & RNase Free	Not Applicable.

BSE/TSE	No use of any raw material produced from or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE).
Cytotoxicity	Not Applicable.
Non-Pyrogenic Statement	Not Applicable.
Latex Statement	Not Applicable.
BPA Statement	Bisphenol-A (BPA) and Bisphenol-F (BPF) are not used in the manufacture or the formulation of the raw materials and are not expected to be present.
DEHP Statement	Not Applicable.
RoHS	No substances (Lead, Cadmium, Mercury, Hexavalent Chromium (Cr ⁶⁺), Poly Brominated Biphenyls (PBB), Poly Brominated Diphenyl ethers (PBDE), Benzyl Butyl Phthalate (BBP), Dibutyl Phthalate (DBP), Diisobutyl Phthalate (DIBP)) are used in manufacturing the raw materials and final product. No routine analysis is performed.
REACH Statement	Not Applicable.
Storage Conditions	Store at room temperature.
Shelf Life	4 years for sterile product and Indefinite shelf life for non-sterile product.

Disclaimer: VWR'S Quality Management System (QMS) is ISO 9001:2015 certified as a wholesale distributor of finished products, ensuring effective oversight of finished product handling, traceability and distribution processes. VWR is not a manufacturer and does not engage in manufacturing, co- manufacturing, or processing activities such as design, packaging, labelling, testing, or quality release. Consequently, the supplier or manufacturer's QMS certification may differ from VWR's.

VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. This product is intended for research use only. It is not for use in diagnostic procedures, therapeutic applications, or any human or clinical use. This document has been produced electronically and is valid without a signature.

Version History

Sl.no	Version number	Description	Month & Year
1	1.0	Created General CoQ using new layout	April 2026