

Certificate of Conformity

VWR Catalogue Number	See List
Description	Flasks, flat bottom, narrow or wide neck
Country of Origin	Manufactured in Czech Republic
Date of Issue (yyyy-mm-dd)	2022-05-17

Narrow neck

Catalogue Number	Capacity [mL]	Height [mm]	Neck Ø ext [mm]	Ø ext [mm]
214-1187	100	110	22	64
214-1188	250	140	34	85
214-1189	500	170	34	105
214-1190	1000	200	42	131
214-1191	2000	250	50	166

Wide neck

Catalogue Number	Capacity [mL]	Height [mm]	Neck Ø ext [mm]	Ø ext [mm]
214-1192	100	110	34	64
214-1193	250	140	50	85
214-1194	500	170	50	105
214-1195	1000	200	50	131
214-1196	2000	250	76	166

Quality System Compliance

Products are manufactured under the **ISO 9001** certified Laboratory. Products are inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable QC requirements.

Product Specifications

Material	ISO 3585 Borosilicate glass 3.3, clear.
Sterilization	Hot air sterilization, in the oven up to 140 °C Steam sterilization, in an autoclave 121 °C/ 20 min 134 °C/ 10 min.
ATP Assay	Not Applicable.
DNase & RNase Free	This product is free of any detectable DNase/RNase contamination.
BSE/TSE	No use of any raw material produced from or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE).
Cytotoxicity	Not Applicable.
Non-Pyrogenic Statement	Not Applicable.
Latex Statement	This product is not made from natural rubber latex.
BPA Statement	Bisphenol are not used in the manufacture of the raw material and are not expected to be present.
DEHP Statement	Not Applicable.

RoHS	No substances (Lead, Cadmium, Mercury, Hexavalent Chromium (Cr6+), Poly Brominated Biphenyls (PBB), Poly Brominated Diphenyl ethers (PBDE), Bis(2-Ethyhexyl) Phthalate (DEHP), Benzyl Butyl Phthalate (BBP), Dibutyl Phthalate (DBP), Diisobutyl Phthalate (DIBP)) are used in manufacturing the raw materials and final product. No routinely analyse is performed.
REACH Statement	No substances as listed on >Substance of Very High Concern< candidate list published 2022-01-17 are used in manufacturing of the raw materials used to produce this product. No routinely analysis is performed.
Storage Conditions	Store at room temperature in normal warehouse conditions; Protect from any possible contamination.
Shelf Life	No predefined shelf life.

Disclaimer: VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. The purpose of the product is for use in laboratory only.
This document has been produced electronically and is valid without a signature.