

Certificate of Conformity

VWR Catalogue Number	See List
Description	Zip-Seal Bags
Country of Origin	Manufactured in China
Date of Issue (yyyy-mm-dd)	2025-07-11

Single Pocket, Non-Sterile

Catalogue Number	Dimension [mm]
129-0293	60 × 40
129-0294	80 × 60
129-0309	100 × 70
129-0295	120 × 80
129-0296	150 × 100
129-0297	170 × 120
129-0298	220 × 160
129-0299	250 × 180
129-0306	300 × 200
129-0307	350 × 250
129-0308	400 × 300

Single Pocket with writing area, Non-Sterile

Catalogue Number	Dimension [mm]
129-0310	60 × 40
129-0311	80 × 60
129-0320	100 × 70
129-0312	120 × 80
129-0313	150 × 100
129-0314	170 × 120
129-0315	220 × 160

Catalogue Number	Dimension [mm]
129-0316	250 × 180
129-0317	300 × 200
129-0318	350 × 250
129-0319	400 × 300

Single Pocket with writing area, Sterile

Catalogue Number	Dimension [mm]
129-0325	300 × 200

Double Pocket, Non-Sterile

Catalogue Number	Dimension [mm]
129-0321	170 × 120
129-0322	220 × 160
129-0323	250 × 180
129-0324	300 × 200

Quality System Compliance

This certificate confirms that the product has been sourced from a supplier or manufacturer who has declared their Quality Management System (QMS) complies with **ISO 9001:2015**.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable QC requirements. Inspected according to AQL level II 4.0 sampling plan per MIL-STD-105E.

Weld strength : $\geq 15\text{N}/25\text{ mm}$.
Tensile Strength : (Vertical $>10\text{MP}/25\text{ mm}$ / Horizontal $> 20\text{MP}/25\text{ mm}$.
Specific Elongation : (Vertical $\geq 450\%$ / Horizontal $\geq 250\%$.

Product Specifications

Material Low Density Polyethylene (LDPE) C150Y + Linear Low-Density Polyethylene (LDPE) 218W (virgin).

Monomers and additives used to produce the products are listed in List of Authorised Substances of Regulation EU/10/2011.

The material grades conform to EU Directive 1935/2004 on materials and articles intended to come into contact with food.

Sterilization	Product labeled as sterile is gamma irradiated and dose released upon ISO11137 recommended practices in effect at the time of validation and meets a minimum requirement of 10 ⁻⁶ SAL Which is not applicable for Non-Sterile products.
ATP Assay	Not applicable.
DNase & RNase Free	Not applicable.
BSE/TSE	No use of any raw material produced from or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE).
Cytotoxicity	Not applicable.
Non-Pyrogenic Statement	Not applicable.
Latex Statement	This product is not made from natural rubber latex.
BPA Statement	Not applicable.
DEHP Statement	Not applicable.
RoHS	Not applicable.
REACH Statement	Not applicable.
Storage Conditions	Room temperature.
Shelf Life	Products labelled as sterile has shelf life period of 3 years which is not applicable for non-sterile products.

Disclaimer: VWR's Quality Management System (QMS) is ISO 9001:2015 certified as a wholesale distributor of finished products, ensuring effective oversight of finished product handling, traceability, and distribution processes. VWR is not a manufacture and does not engage in manufacturing, co-manufacturing or processing activities such as design, packaging, labeling, testing or quality release. Consequently, the supplier or manufacturer's QMS certification may differ from VWR's.

VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. This product is intended for research use only. It is not for use in diagnostic procedures, therapeutic applications, or any human or clinical use.

This document has been produced electronically and is valid without a signature.

Version History

Sl.no	Version number	Description	Month & Year
1	1.0	Created General CoQ using new layout	08-2021
3	2.0	Modified Quality system compliance and disclaimer	07-2025