




The Multiroir range of presentation trays is designed to meet the needs of laboratories. These trays combine elegant design with functionality.

Made from soft-touch polypropylene, they are available in six colors: purple, orange, blue, translucent, pink, and green. Their design ensures high impact resistance and optimal stacking, with external dimensions of 240 x 167 x H 52 mm and internal dimensions of 230 x 158 x H 48 mm, perfectly suited for tracking sheets and other laboratory documents.

PRODUCT ADVANTAGES

- **High-Quality material** : Made from polypropylene with a soft touch, they combine robustness and aesthetic appeal.
- **Sleek design** : their design is crafted to meet the demands of optical professionals, blending elegance with practicality.
- **Shock resistance** : their high shock resistance ensures increased durability and longevity.
- **Efficient stacking**: designed for optimal stacking, they provide space-saving benefits and easy organization.

SPECIFICATIONS

SKU	MULC26506PENPR
Ext. dimensions	240 x 167 x 52 mm
Int. dimensions	230 x 158 x 48 mm
Weight	180 g
Material	Polypropylene
Wall type / bottom	Solid
Color	Pink
Recyclable	Yes
Food use	Yes
Shipping lead time	5 working days
Origin	France