

Certificate of Quality

VWR Catalogue Number	58948-072
Description	VWR® Spinbar® Polygon Stir Bars, Assorted, PTFE Coating, Pack of 18
Country of Origin	Manufactured in USA
Date of Issue (yyyy-mm-dd)	2025-07-08

Quality System Compliance

This certificate confirms that the product has been sourced from a supplier or manufacturer who has declared their Quality Management System (QMS) complies with [ISO 9001:2015](#) standard.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable QC requirements.

Autoclavable @121° C; for 20 minutes
Tolerance +/- 0.020"

Product Specifications

Material	Dupont White Teflon® #807N-X, ALNICO 5 Magnet, Polystyrene Storage Box Materials used meet FDA Regulations 21CFR177.1550, 21CFR177.1520 and Teflon used to totally encapsulate the magnetic bar meets US Pharmacopeias Class VI. This product does not contain or come in contact Formaldehyde or other contaminants. Includes stir bars in seven sizes: 20 x 7mm ($25/32 \times 9/32$"), 25 x 8mm ($31/32 \times 5/16$"), 30 x 8mm ($13/16 \times 5/16$"), 40 x 8mm ($19/16 \times 5/16$"), 50 x 8mm ($131/32 \times 5/16$"), 60 x 8mm ($23/8 \times 5/16$"), and 80 x 10mm ($35/32 \times 13/32$")
Sterilization	The product is non-sterile.
ATP Assay	Not applicable.
DNase & RNase Free	Not applicable.

BSE/TSE	No use of any raw material produced from or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE).
Cytotoxicity	Not applicable
Non-Pyrogenic Statement	Not applicable
Latex Statement	The product is free from Latex.
BPA Statement	Bisphenols are not used in the manufacture of the raw material and are not expected to be present.
DEHP Statement	Diethylhexyl phthalate (DEHP) is not included in the manufacture or the formulation of the raw material and are not expected to be present.
RoHS	No substances (Lead, Cadmium, Mercury, Hexavalent Chromium (Cr ⁶⁺), Poly Brominated Biphenyls (PBB), Poly Brominated Diphenyl ethers (PBDE), Benzyl Butyl Phthalate (BBP), Dibutyl Phthalate (DBP), Diisobutyl Phthalate (DIBP)) are used in manufacturing the raw materials and final product. No routinely analyse is performed.
REACH Statement	No substances as listed on >Substance of Very High Concern< candidate list are used in manufacturing of the raw materials used to produce this product. No routinely analyse is performed.
Storage Conditions	Store at low humidity, out of direct sunlight at ambient warehouse temperature, in low humidity conditions, protected from direct sunlight and away from harsh chemicals.
Shelf Life	There is no definite shelf life or expiration for the product. The product exhibits no degradation for 5 years if maintained in above storage condition.

Disclaimer: VWR'S Quality Management System (QMS) is ISO 9001:2015 certified as a wholesale distributor of finished products, ensuring effective oversight of finished product handling, traceability and distribution processes. VWR is not a manufacturer and does not engage in manufacturing, co- manufacturing, or processing activities such as design, packaging, labelling, testing, or quality release. Consequently, the supplier or manufacturer's QMS certification may differ from VWR's.

VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. This product is intended for research use only. It is not for use in diagnostic procedures, therapeutic applications, or any human or clinical use. This document has been produced electronically and is valid without a signature.

Version History

Sl.no	Version number	Description	Month & Year
1	1.0	Created General CoQ using new layout	July 2025