

Material No.: GT20-05
Revision No.: 0

Product Specification

GMP Manufactured Product

Test	Specification
Appearance (Clear, colorless)	Passes Test
Identification (SDS Page)	Passes Test
Purity by SDS-PAGE	$\geq 99 \%$
Activity (Units/uL)	≥ 250
Specific activity (Units/mg)	≥ 1000000
Endotoxin (EU/1000U)	≤ 0.25
TAMC (cfu/100,000U)	≤ 10
TYMC (cfu/100,000U)	≤ 10
Protease Activity	None Detected
Mycoplasma	None Detected
Adventitious viruses	None Detected

GMP Manufactured Product
For drug manufacturing or R&D use only.
Storage Condition: Store at or below -20°C

Packaging Site: Phillipsburg Mfg Ctr & DC