

Certificate of Quality

VWR Catalogue Number	See List
Description	VWR® Disposable Scrub Shirts with long sleeve, adjustable neck
Country of Origin	Manufactured in USA
Date of Issue (yyyy-mm-dd)	2021-08-16

VWR Catalogue Number	Colour	Clothing size
10755-052	Blue	S
10754-966	Blue	M
10754-968	Blue	L
10754-970	Blue	XL
10754-972	Blue	XXL
10754-974	Blue	3XL
10754-976	Blue	4XL
10754-978	Blue	5XL

Quality System Compliance

Products are manufactured in accordance with the manufacture's quality system. Products are inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable product specifications. The product meets Personal Protective Equipment (PPE) requirements.

Product Specifications

Material	Spunbond–Meltblown–Spunbond (SMS) Fabric
Sterilization	Not applicable.
ATP Assay	Not applicable.
DNase & RNase Free	This product is free of any detectable DNase/RNase contamination.

BSE/TSE	No use of any raw material produced from or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE).
Cytotoxicity	Testing is conducted to quality all material resins and have been shown to be non-toxic.
Non-Pyrogenic Statement	The product is free from pyrogens and endotoxins.
Latex Statement	The product is Latex free.
BPA Statement	Bisphenol are not used in the manufacture of the raw material and are not expected to be present.
DEHP Statement	Not Applicable.
RoHS	Not Applicable.
REACH Statement	Not Applicable.
Storage Conditions	Store at Ambient condition.
Shelf Life	10 years.

Disclaimer: VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. The purpose of the product is for use in laboratory only.
This document has been produced electronically and is valid without a signature.