

Product Specification Sheet

VWR DISPOSABLE STRAW PIPETTE, STERILE, WITHOUT GRADUATION

VWR International certifies that the specified products meet all VWR specifications and quality criteria for release. The products comply in every respect with the relevant particular specifications, drawings, relevant standards and regulations in force.

VWR International certifie que les produits mentionnés respectent toutes les spécifications et critères de qualité VWR. Les produits satisfont aux spécifications particulières, aux plans, aux normes et règlements en vigueur s'y rapportant.

VWR International bestätigt hiermit, daß die angegebenen Produkte alle VWR Spezifikationen und Qualitätskriterien zur Freigabe erfüllen. Die Produkte entsprechen in jeder Hinsicht den maßgeblichen Spezifikationen, den Zeichnungen sowie den relevanten Normen und gültigen Vorschriften.

Cat. No	Nominal capacity	Length	Packaging
612-1675	1 ml	190 mm	40 bags of 25 pieces

Material: **Polypropylene**

Storage Condition (°C): **Room temperature (from 15 to 30°C)**

Sterility Assurance Level based on ISO 11137: **Sterile SAL 10⁻³**

Physical properties: **Temperature resistance main body -20°C / + 120°C**

Flexible to rigid, good mechanical resistance, high chemical resistance

Other Characteristics:

- **Traceability including the expiry date on the secondary packaging**
- **Zip bag can be pressed closed thus ensuring the protection of the products**
- **Tamper-evident seal guarantees the sterility of the product when it is first opened**

Disclaimer: VWR wants to state that these specifications can only be regarded as general information and will not discharge the user from his obligation to ensure the product is suitable for the intended use - If there is any doubt ask us.

Date of issue: **May 2013**