

Material	20808.291
Material description	Butan-1-ol
Grade	GPR RECTAPUR
CAS Number	71-36-3
Molecular formula	$\text{H}_3\text{C}(\text{CH}_2)_2\text{CH}_2\text{OH}$
Molecular mass	74.12

Characteristics	Specifications
Assay	$\geq 98.5 \%$
Boiling point	116.6 - 118.3 °C
Density (20/4)	0.807 - 0.812
Evaporation residue	$\leq 100 \text{ ppm}$

Signature

We certify that this batch conforms to the specifications listed above.

This document has been produced electronically and is valid without a signature.

Isabelle Habay, Head of Laboratory - Briare
VWR International S.A.S.; Z.I. de Vaugereau; FR-45250
Briare; France

For Professional use in Laboratory or Manufacturing. Not for use as an Active Pharmaceutical Ingredient or Food or Animal Feed. Suitability and intended use of the product remains the responsibility of the user.