


Recombinant E. coli RuvA Protein (Functional) (A458)

Specifications:

Name:	Recombinant E. coli RuvA Protein (Functional)
Applications:	SDS-PAGE, Functional Studies
Applications:	Functional as Holliday junction specific binding protein, which promotes Holliday-junction branch migration in combination with RuvB protein.
Expression System:	Escherichia coli
Nature:	Recombinant
Protein Species:	Escherichia coli
Protein Length:	Full length protein
Molecular Weight:	22 kDa
Biologically Active:	Yes
Purity:	> 90% (by SDS-PAGE).
Product Form:	Liquid
Concentration:	1.0-2.7 mg/ml
Formulation:	Supplied in 10 mM Tris-HCl, pH 7.5, with 2 mM EDTA, 100 mM NaCl, 5 mM Mercaptoethanol, and 50% Glycerol.
Storage:	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
Disclaimer:	This product is for research use only. It is not intended for diagnostic or therapeutic use.

Images:



Recombinant E. coli RuvA Protein (Functional)