

Material	28726.320
Material description	Tributyl phosphate
Grade	GPR RECTAPUR

CAS Number	126-73-8
Molecular formula	$(\text{CH}_3(\text{CH}_2)_3\text{O})_3\text{PO}$
Molecular mass	266.32

Characteristics	Specifications
Assay	≥ 99.0 %
Appearance	Clear liquid
Colouration	≤ 50 APHA
Density (20/4)	0.975 - 0.980
Refractive index (20°C)	1.423 - 1.425
Water	≤ 0.2 %

Signature

We certify that this batch conforms to the specifications listed above.

This document has been produced electronically and is valid without a signature.

Isabelle Habay, Head of Laboratory - Briare
VWR International S.A.S.; Z.I. de Vaugereau; FR-45250
Briare; France

For Professional use in Laboratory or Manufacturing. Not for use as an Active Pharmaceutical Ingredient or Food or Animal Feed. Suitability and intended use of the product remains the responsibility of the user.