



Specification

1.06462.1000 Sodium hydroxide pellets EMPLURA®

Specification		
Assay (acidimetric, NaOH)	≥ 97.0	%
Carbonate (as Na ₂ CO ₃)	≤ 1.0	%
Chloride (Cl)	≤ 0.012	%
Sulfate (SO ₄)	≤ 0.010	%
Heavy metals (as Pb)	≤ 0.002	%
Al (Aluminium)	≤ 0.002	%
Fe (Iron)	≤ 0.002	%

Dr. Dimitrij Ryvlin

Responsible laboratory manager quality control

This document has been produced electronically and is valid without a signature.