

Technical Data Sheet

VWR Catalogue Number	See List
Description	VWR® High Performance transfer pipet
Country of Origin	Manufactured in Italie
Date of Issue (yyyy-mm-dd)	2021-04-19

GRADUATED



Catalogue Number	Capacity ml	Drops per ml	Bulb draw ml	Stem Ø mm	Length mm	Grad up to ml	Sterile	Packed	PK
612-1681	7	21	3,5	7.8	150	3	-	Bulk	500
612-1682	7	21	3,5	7.8	150	3	+	5 per bag	500
612-1683	7	21	3,5	7.8	150	3	+	10 per bag	500
612-1747	7	21	3,5	7.8	150	3	+	Individual	500



Catalogue Number	Capacity ml	Drops per ml	Bulb draw ml	Stem Ø mm	Length mm	Grad up to ml	Sterile	Packed	PK
612-1684	5	23	3,5	5	150	1	-	Bulk	500
612-1685	5	23	3,5	5	150	1	+	Individual	500
612-1686	5	23	3,5	5	150	1	+	5 per bag	500
612-1687	5	23	3,5	5	150	1	+	10 per bag	500

Technical Data Sheet

shown at 80% of actual size



Catalogue Number	Capacity ml	Drops per ml	Bulb draw ml	Stem Ø mm	Length mm	Grad up to ml	Sterile	Packed	PK
612-2842	6	22	3,5	5	230	1.5	-	Bulk	400
612-2843	6	22	3,5	5	230	1.5	+	Individual	500
612-2844	6	22	3,5	5	230	1.5	+	5 per bag	500
612-2845	6	22	3,5	5	230	1.5	+	10 per bag	500



Catalogue Number	Capacity ml	Drops per ml	Bulb draw ml	Stem Ø mm	Length mm	Grad up to ml	Sterile	Packed	PK
612-2851	3	22	1,5	5	138	1	-	Bulk	500

Technical Data Sheet

NARROW STEM



Catalogue Number	Capacity (ml)	Drops per ml	Bulb draw ml	Stem Ø (mm)	Length (mm)	Sterile	Packed	Pk
612-1688	4	28	3.5	2.5	150	-	Bulk	500
612-1756	4	28	3.5	2.5	150	+	Individual	500
612-1689	4	28	3.5	2.5	150	+	5 per bag	500
612-1690	4	28	3.5	2.5	150	+	10 per bag	500



Catalogue Number	Capacity (ml)	Drops per ml	Bulb draw ml	Stem Ø (mm)	Length (mm)	Sterile	Packed	Pk
612-2841	4	22	3.5	2.5	83	-	Bulk	1.000



Catalogue Number	Capacity (ml)	Drops per ml	Bulb draw ml	Stem Ø (mm)	Length (mm)	Sterile	Packed	Pk
612-2848	1.2	25	0.9	2.5	63	-	Bulk	500

BLOOD BANK



Catalogue Number	Capacity (ml)	Drops per ml	Bulb draw ml	Stem Ø (mm)	Length (mm)	Sterile	Packed	Pk
612-2850	5	22	2,4	6,7	155	-	Bulk	500

Technical Data Sheet

FINE TIP PIPET



Catalogue Number	Capacity ml	Drops per ml	Bulb draw ml	Stem Ø mm	Length mm	Sterile	Packed	PK
612-2853	1,5	50	1	3	104	-	Bulk	400



Catalogue Number	Capacity ml	Drops per ml	Bulb draw ml	Stem Ø mm	Length mm	Sterile	Packed	PK
612-2854	7	66	6	6,3	150	+	/ 5	500
612-2855	7	66	6	6,3	150	+	/ 10	500
612-4464	7	66	6	6,3	150	-	Bulk	400
612-4465	7	66	6	6,3	150	+	/ 1	500



Catalogue Number	Capacity ml	Drops per ml	Bulb draw ml	Stem Ø mm	Length mm	Sterile	Packed	PK
612-2856	5	65	3,5	5	144	-	Bulk	500
612-2857	5	65	3,5	5	144	+	/ 1	500
612-2858	5	65	3,5	5	144	+	/ 10	500

Technical Data Sheet

PADDLE PIPET




Catalogue Number	Capacity ml	Drops per ml	Bulb draw ml	Stem Ø mm	Length mm	Sterile	Packed	PK
612-2852	3	33	0,2	3,8	130	-	Bulk	500

EXTRA LARGE



Catalogue Number	Capacity	Dops/ml	Bulb draw ml	Stem Ø mm	Length mm	Packed	Pack Size
612-2847	14.0 ml	18	10	4.0	170	Bulk	200



Catalogue Number	Capacity	Dops/ml	Bulb draw ml	Stem Ø mm	Length mm	Packed	Pack Size
612-2859	23.0 ml	22	10.0	9.0	300	Bulk	600

Quality System Compliance

Products are inspected and controlled in accordance with the current applicable product specification. Inspection records are reviewed and maintained.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable QC requirements

Product Specifications

Material LDPE
Packaging : Paper-film pouch

Technical Data Sheet

Top Sheet Material : Medical Grade Paper Plus 60 GR
Bottom Sheet Material : PA/PE Film - 100µ

Sterilization	Sterilisation method: have been EO treated (ISO 11135)
Non-Pyrogenic Statement	Not tested
ATP Assay	Not tested
DNase & RNase Free	Not tested
BSE/TSE	Resin suppliers have indicated that the LDPE resin is produced with chemicals that are synthetically derived and is free from animal derived materials, including BSE/TSE (Bovine Spongiform Encephalopathy/Transmissible spongiform Encephalopathy).
Cytotoxicity	Not tested
Latex Statement	Latex-free
BPA Statement	Not tested
DEHP Statement	Not tested
RoHS	The resin suppliers have indicated that the material they supply is compliant with revised RoHS directive 2011/65/EC and RoHS 2 restricting the use of heavy metals, PBB'S, and PBDE's. This statement does not constitute legal or business advice. This is provided for informational purposes only; and is based on information provided by our suppliers at this time this document was published. This compliance information is subject to change at any time.
REACH Statement	The resin suppliers have indicated that the material they supply is compliant with the REACH Directives 552/2009, 1907/2006, 197 SVHC and do not contain any SVHCs in amounts greater than 0.1% (w/w). This statement does not constitute legal or business advice. This is provided for informational purposes only; and is based on information provided by our suppliers at this time this document was published. This compliance information is subject to change at any time.
Storage Conditions	No special storage conditions required
Shelf Life	Not applicable

Disclaimer: VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. The purpose of the product is for use in laboratory only.

This document has been produced electronically and is valid without a signature