

Calibration solutions for conductivity – WTW – WTWA300580



- [Specifications](#)

Solutions for the calibration of standard and pure water measuring cells

Calibration solutions for checking and calibrating of a system consisting of a meter and a conductivity cell

- Optimal reliability by PTB/NIST traceable standard solutions
- Special calibration solutions for pure water cells
- High accuracy

E-Set Trace: PTB/NIST–traceable conductivity standard, 0,01 mol/l KCl, 1413 µS/cm at 25 °C, accuracy +/- 0.5 %, 6 bottles with 50 ml each

KS 100 µS: NIST traceable KCl containing standard solution, accuracy +/- 3 %, bottle with 250 ml

KS 5 µS: NIST traceable KCl containing standard solution, accuracy +/- 2 %, bottle with 250 ml