



Specification

1.02239.0500 Peptone from casein pancreatically digested free from sulfonamide antagonists for microbiology

Specification

Appearance	Lightly yellow to lightly brownishyellow fine powder. Characteristic, not rotten smell.	
Identity	passes test	
pH-value (5 %; water)	5.0 - 6.0	
Total nitrogen (N) (Kjeldahl)	13.8 - 15.8	%
Amino nitrogen (as N)	4.7 - 6.7	%
Sulfated ash (800 °C)	≤ 5.0	%
Loss on Drying (105 °C)	≤ 6.5	%
Suitability for microbiology	passes test	

Stefanie Fischer

Responsible laboratory manager quality control

This document has been produced electronically and is valid without a signature.