

Certificate of Conformity

VWR Catalogue Number	113-9074
Description	Sleeve protectors, VWR® Basic, 400 mm Length, 20 µm Thickness, Blue
Lot Number	20211201
Date of Manufacture (yyyy-mm-dd)	2021-12-01
Country of Origin	Manufactured in China
Date of Issue (yyyy-mm-dd)	2023-02-06

Quality System Compliance

Products are manufactured under the **ISO 9001:2015 & ISO 13485:2016** standard. Products are inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP.

QC Testing

Representative products samples are collected and inspected according to 4.0 AQL level II sampling plan

	Test / Standard	MD	CD
Break Tensile Properties [N]	ISO 13934-1	23.41	11.52
Elongation Rate [%]	ISO 13934-1	195.18	142.33

Product Specifications

Material	PE Film.
Sterilization	Non-Sterile.
ATP Assay	Not Applicable.
DNase & RNase Free	Not Applicable.

BSE/TSE	No use of any raw material produced from or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE).
Cytotoxicity	Not Applicable.
Non-Pyrogenic Statement	Not Applicable.
Latex Statement	Nature rubber latex but without exposure to outside.
BPA Statement	Not Applicable.
DEHP Statement	Not Applicable.
RoHS	Not Applicable.
REACH Statement	Not Applicable.
Storage Conditions	Temp: 0- 50°C, Humidity: 45-85% without sunshine exposure.
Shelf Life	Not Applicable.

Disclaimer: VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. The purpose of the product is for use in laboratory only.
This document has been produced electronically and is valid without a signature.

Version History

Sl.no	Version number	Description	Month & Year
1	1.0	Created General CoQ using new layout	02-2023