

Certificate of Conformity

VWR Catalogue Number	See List
Description	Disposable Funnels with Short Stem
Country of Origin	Manufactured in India
Date of Issue (yyyy-mm-dd)	2021-11-24

Polypropylene

Catalogue Number	Funnel Outer Diameter [mm]	Length [mm]	Stem Outer Diameter [mm]	Stem length [mm]
221-0251	65	74.4	7.4/4.3	22.7
221-0252	55	71.3	7.4/4.3	27.1

Polystyrene

Catalogue Number	Funnel Outer Diameter [mm]	Length [mm]	Stem Outer Diameter [mm]	Stem length [mm]
221-0253	55	71.3	7.4/4.3	27.1

Quality System Compliance

Products are manufactured under the **ISO 9001:2015** and **ISO 13485:2016** standards. Products are inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable QC requirements.

For 221-0251 and 221-0252

Test Name	Test Reference	Accuracy	Result
Autoclaving	QWI-QC-07	General Insp. Level-1, AQL -1.5	PASS
Dimensional	QSF-QCIR-06	General Insp. Level-1, AQL -1.5	Within Tolerance

For 221-0253

Test Name	Test Reference	Accuracy	Result
Autoclaving	QWI-QC-07	General Insp. Level-1, AQL -1.5	NA
Dimensional	QSF-QCIR-06	General Insp. Level-1, AQL -1.5	Within Tolerance

Product Specifications

Material	Polypropylene (PP), 100 % Virgin, FDA 21CFR Compliance. or Polystyrene (PS) 100 % Virgin, FDA 21CFR Compliance.
Sterilization	Not Applicable.
ATP Assay	Not Applicable.
DNase & RNase Free	Not Applicable.
BSE/TSE	No use of any raw material produced from or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE).
Cytotoxicity	Not Applicable.
Non-Pyrogenic Statement	Not Applicable.
Latex Statement	This product is not made from natural rubber latex.
BPA Statement	Bisphenol are not used in the manufacture of the raw material and are not expected to be present..
DEHP Statement	Please see RoHS information below.
RoHS	No substances (Lead, Cadmium, Mercury, Hexavalent Chromium (Cr6+), Poly Brominated Biphenyls (PBB), Poly Brominated Diphenyl ethers (PBDE), Bis(2-Ethyhexyl) Phthalate (DEHP), Benzyl Butyl Phthalate (BBP), Dibutyl Phthalate (DBP), Diisobutyl Phthalate (DIBP)) are used in manufacturing the raw materials and final product. No routinely analyse is performed.

REACH Statement	No substances as listed on >Substance of Very High Concern< candidate list published 2022-01-17 are used in manufacturing of the raw materials used to produce this product. No routinely analysis is performed.
Storage Conditions	Store at Ambient temperature.
Shelf Life	Not Applicable.

Disclaimer: VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. The purpose of the product is for use in laboratory only.
This document has been produced electronically and is valid without a signature.