

# Certificate of Conformity

|                                         |                                                   |
|-----------------------------------------|---------------------------------------------------|
| <b>VWR Catalogue Number</b>             | VWRA734-3330                                      |
| <b>Description</b>                      | VWR®, Classic, PCR 96-well plate sealing membrane |
| <b>Lot Number</b>                       | 20707/2526_12                                     |
| <b>Date of Manufacture (yyyy-mm-dd)</b> | 2025-12                                           |
| <b>Country of Origin</b>                | Manufactured in India                             |
| <b>Date of Issue (yyyy-mm-dd)</b>       | 2026-05-05                                        |

## Quality System Compliance

This certificate confirms that the product has been sourced from a supplier or manufacturer who has declared their Quality Management System (QMS) complies with [ISO 9001:2015 & ISO ISO 13485:2016](#).

## QC Testing

Representative products samples are collected and inspected in accordance with current applicable QC requirements.

## Product Specifications

|                               |                                                                                                                                                                            |
|-------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Material</b>               | Polypropylene (PP).<br><br>Resins are USP Class-VI certified, FDA regulation CFR 21, and are FDA approved for food contact.                                                |
| <b>Sterilization</b>          | Not Applicable.                                                                                                                                                            |
| <b>ATP Assay</b>              | Not Applicable.                                                                                                                                                            |
| <b>DNase &amp; RNase Free</b> | This product is free of any detectable DNase/RNase contamination. Test limit is $1.0 \times 10^{-6}$ Kunitz Units for DNase & $1.0 \times 10^{-9}$ Kunitz Units for RNase. |

|                                |                                                                                                                                                                                                                                                                                                                                                                                                                                  |
|--------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>BSE/TSE</b>                 | No use of any raw material produced from or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE) as defined in EC Directive 1069/2009/EC and EC Directive 142/2011/EC. |
| <b>Cytotoxicity</b>            | This product is free of cytotoxicity effects.                                                                                                                                                                                                                                                                                                                                                                                    |
| <b>Non-Pyrogenic Statement</b> | Endotoxin levels are tested by Limulus Amoebocyte Lysate (LAL) Assay. Product is free from Endotoxin. The acceptance level for product is $\leq 0.01$ EU/ ml (LAL Gel Clot Method).                                                                                                                                                                                                                                              |
| <b>Latex Statement</b>         | Product is free from latex.                                                                                                                                                                                                                                                                                                                                                                                                      |
| <b>BPA Statement</b>           | Bisphenols has not been intentionally incorporated in the manufacturing of the products.                                                                                                                                                                                                                                                                                                                                         |
| <b>DEHP Statement</b>          | DEHP has not been intentionally incorporated in the manufacturing of the products.                                                                                                                                                                                                                                                                                                                                               |
| <b>RoHS</b>                    | No substances (Lead, Cadmium, Mercury, Hexavalent Chromium ( $\text{Cr}^{6+}$ ), Poly Brominated Biphenyls (PBB), Poly Brominated Diphenyl ethers (PBDE), Benzyl Butyl Phthalate (BBP), Bis(2-Ethylhexyl) Phthalate (DEHP), Dibutyl Phthalate (DBP), Diisobutyl Phthalate (DIBP)) are used in manufacturing the raw materials and final product. No routinely analyse is performed.                                              |
| <b>REACH Statement</b>         | No substances as listed on >Substance of Very High Concern< candidate list are used in manufacturing of the raw materials used to produce this product. No routinely analysis is performed.                                                                                                                                                                                                                                      |
| <b>Storage Conditions</b>      | Ambient Room temperature.                                                                                                                                                                                                                                                                                                                                                                                                        |
| <b>Shelf Life</b>              | 3 years.                                                                                                                                                                                                                                                                                                                                                                                                                         |

Disclaimer: VWR's Quality Management System (QMS) is ISO 9001:2015 certified as a wholesale distributor of finished products, ensuring effective oversight of finished product handling, traceability, and distribution processes. VWR is not a manufacture and does not engage in manufacturing, co-manufacturing or processing activities such as design, packaging, labeling, testing or quality release. Consequently, the supplier or manufacturer's QMS certification may differ from VWR's.

VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. This product is intended for research use only. It is not for use in diagnostic procedures, therapeutic applications, or any human or clinical use.

This document has been produced electronically and is valid without a signature.