

Product Specification

Meets ACS Specifications.

Test	Specification
Assay	$\geq 99.9 \%$
Clarity of Solution	passes test
Free Acid	passes test
Lead (Pb)	$\leq 0.001 \%$
Substances not Precipitated by HCl	$\leq 0.01 \%$
Sulphate (SO ₄)	$\leq 0.002 \%$
Trace Impurities (in ppm)	
Chloride (Cl)	≤ 5
Copper (Cu)	≤ 2
Iron (Fe)	≤ 2

For Laboratory, Research or Manufacturing use.



Barbara Rykowska
Quality Assurance Coordinator

For questions on this Product Specification please contact Technical Service at technical.service@avantorsciences.com
Avantor Performance Materials Poland S.A. Sowinskiego 11 Gliwice 44-101 Poland Tel: +48 322392000