

Certificate of Quality

VWR Catalogue Number	See List
Description	VWR® Sharps Containers with Rotating Lid
Country of Origin	Manufactured in USA
Date of Issue (yyyy-mm-dd)	2022-04-08

Rotor Lids Openings

VWR Catalogue Number	Length x Width x Height [inches]	Capacity [gal]	Colour
76319-806	10.20 × 7.0 × 5.60	1	Red
76319-808	10.20 × 7.0 × 10.10	2	Red
76319-804	6.25 × 4.5 × 7.00	0.457	Red

Rotor + Half Hinged Openings

VWR Catalogue Number	Length x Width x Height [inches]	Capacity [gal]	Colour
76319-810	10.30 × 7.0 × 10.10	3	Red
76319-858	17.75 × 15.5 × 17.75	8	Red

Chimney Top Openings

VWR Catalogue Number	Length x Width x Height [inches]	Capacity [gal]	Colour
76319-816	10.20 × 7.0 × 7.00	1	Red
76319-818	10.20 × 7.0 × 12.00	2	Red

Quality System Compliance

Products are manufactured under the **ISO 9001:2015 & ISO 13485:2003** standard. Products are inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable product specifications.

Dimensional Analysis: ASQ Z1.4-2003, General Insp. Level- II, AQL- 0.65. **Result:** Within tolerance

Product Specifications

This product is Class I category device as defined by the FDA in 21CFR Parts 862-892.

Material	Container: High-Density Polyethylene (HDPE), 100% Virgin, FDA 21 CFR Compliance, USP Class VI passed Lid: Blend of High-Density Polyethylene (HDPE) and Polyolefin Elastomer, 100% Virgin, FDA 21 CFR Compliance, USP Class VI passed
Sterilization	Not Applicable.
ATP Assay	Not Applicable.
DNase & RNase Free	Not Applicable.
BSE/TSE	No use of any raw material produced from, or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE).
Cytotoxicity	Not Applicable.
Non-Pyrogenic Statement	Not Applicable
Latex Statement	Not Applicable.
BPA Statement	Not Applicable.
DEHP Statement	Not Applicable.
RoHS	Not Applicable.
REACH Statement	Not Applicable.
Storage Conditions	Climate controlled.
Shelf Life	Not Applicable.

Disclaimer: VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. The purpose of the product is for use in laboratory only.
This document has been produced electronically and is valid without a signature.