



VWR® 1.1mL Milk testing Pipet

Description: 1.1mL Milk testing Pipet


Purpose: Use for liquid transfer

Materials

Pipet: GPPS (General polystyrene), Compliance with USP Class VI

Filter: PO (Polyolefin)

Color: Clear



Unit: mm

	Dimensions	Tolerance
A	ø 4.5	±0.1
B	ø 2.45	±0.35
C	268.5	±2.0



Features

- Temperature: -20°C~+50°C
- Sterilized by gamma irradiation SAL10⁻⁶(ISO11137)
- Latex Statement: this product is latex free
- DNase/RNase free, non-pyrogenic, non-cytotoxic
- Meet the requirement of BSE/TSE
- Meet the requirement of APHA
- Use dispenser box for easy opening
- Available in Individually packed in peel-to-open wrap
- Autoclavable: No
- Total volume = 1.1 ml
- Graduations are calibrated for accurate dispensing to within ±2%
- Color coded: purple
- Every product is printed with a batch number for traceability
- Shelf Life: 3 years after month of production
- Manufactured in a class 100,000 room environment
- Manufactured under ISO13485:2016 and ISO9001:2015 quality management system

VWR® 1.1mL Milk testing Pipet

VWR EU Cat. No.	Description	Sterile	Qty/Tray	Qty/Pk	Qty/Case
612-5851	1.1mL, Sterilized, Purple, Individually Wrapped, Paper/Plastic	Yes	1	500	500
612-5852	1.1mL, Sterilized, Purple, Plastic pack	Yes		25	1000



VWR® 2.2mL Milk testing Pipet

Description: 2.2mL Milk testing Pipet


Purpose: Use for liquid transfer

Materials

Pipet: GPPS (General polystyrene), Compliance with USP Class VI

Filter: PO (Polyolefin)

Color: Clear



Unit: mm

	Dimensions	Tolerance
A	ø6.3	±0.1
B	ø2.93	±0.35
C	272.0	±2.0



Features

- Temperature: -20°C~+50°C
- Sterilized by gamma irradiation SAL10⁻⁶(ISO11137)
- Latex Statement: this product is latex free
- DNase/RNase free, non-pyrogenic, non-cytotoxic
- Meet the requirement of BSE/TSE
- Meet the requirement of APHA
- Use dispenser box for easy opening
- Available in Individually packed in peel-to-open wrap
- Autoclavable: No
- Total volume = 2.2 ml
- Graduations are calibrated for accurate dispensing to within ±2%
- Color coded: black
- Every product is printed with a batch number for traceability
- Shelf Life: 3 years after month of production
- Manufactured in a class 100,000 room environment
- Manufactured under ISO13485:2016 and ISO9001:2015 quality management system

VWR® 2.2mL Milk testing Pipet

VWR EU Cat. No.	Description	Sterile	Qty/Tray	Qty/Pk	Qty/Case
612-5853	2.2mL, Sterilized, Black, Individually Wrapped, Paper/Plastic	Yes	1	500	500
612-5854	2.2mL, Sterilized, Black, Plastic pack	Yes		25	1000