



VWR® 35mm culture dish with 15mm glass bottom


Description: 35mm VWR culture dish with 15mm glass bottom, sterile

Purpose: Suitable for fluorescent microscope experiment, confocal microscopy and phase-contrast microscopic experiment, etc

Materials

Dish: GPPS (General polystyrene) Color: Clear
Lid: GPPS (General polystyrene) Color: Clear
Glass bottom: D263 M Borosilicate glass Color: Clear

3:Dimensions (Unit: mm)



Features

- Glass bottom diameter of 15mm
- Glass thickness: 0.16mm
- Borosilicate glasses is with high homogeneity, low bubble and inclusion content
- Medical grade shadowless glue, highly transparent achromatic
- Suitable for fluorescent microscope experiment , confocal microscopy and phase-contrast microscopic experiment , etc
- Suitable for all living cells' examinations
- Dish surface is smooth and free from striation to maximize usable area for growth
- The rim on upper side of the lid mates with the dish brim for easy and secure stack
- Flat bottom uniform wall thickness ensures distortion-free bottom
- Lids with several little chimbs to shape vents are available for very effective gas exchange
- Temperature range: -20°C to +50°C
- Autoclavability: No
- Non-pyrogenic
- Sterilisation: sterilized by gamma irradiation SAL10⁻⁶(ISO11137)
- Shelf life: 3 years after month of production



VWR® 35mm culture dish with 15mm glass bottom

VWR NA Cat. No.	VWR EU Cat. No.	Glass Diameter(mm)	Surface Type	Sterile	Qty. per bag/case
75856-740	734-2903	15	Standard, Surface treated	Y	10/240
75856-744	734-2905	15	General, Non-treated	Y	10/240



VWR® 35mm culture dish with 20mm glass bottom

Description: 35mm VWR culture dish with 20mm glass bottom, sterile

Purpose: Suitable for fluorescent microscope experiment, confocal microscopy and phase-contrast microscopic experiment, etc


Materials

Dish: GPPS (General polystyrene) Color: Clear
Lid: GPPS (General polystyrene) Color: Clear
Glass bottom: D263 M Borosilicate glass Color: Clear

Features

- Glass bottom diameter of 20mm
- Glass thickness: 0.16mm
- Borosilicate glasses is with high homogeneity, low bubble and inclusion content
- Medical grade shadowless glue, highly transparent achromatic
- Suitable for fluorescent microscope experiment , confocal microscopy and phase-contrast microscopic experiment , etc
- Suitable for all living cells' examinations
- Dish surface is smooth and free from striation to maximize usable area for growth
- The rim on upper side of the lid mates with the dish brim for easy and secure stack
- Flat bottom uniform wall thickness ensures distortion-free bottom
- Lids with several little chimbs to shape vents are available for very effective gas exchange
- Temperature range: -20°C to +50°C
- Autoclavability: No
- Non-pyrogenic
- Sterilisation: sterilized by gamma irradiation SAL10⁻⁶(ISO11137)
- Shelf life: 3 years after month of production

3:Dimensions (Unit: mm)



VWR® 35mm culture dish with 20mm glass bottom

VWR NA Cat. No.	VWR EU Cat. No.	Glass Diameter(mm)	Surface Type	Sterile	Qty. per bag/case
75856-742	734-2904	20	Standard, Surface treated	Y	10/240
75856-746	734-2906	20	General, Non-treated	Y	10/240