

Certificate of Conformity

VWR Catalogue Number	See List
Description	VWR®, Single-channel Pipettes, Mechanical, Variable Volume
Country of Origin	Manufactured in Poland
Date of Issue (yyyy-mm-dd)	2023-09-20

Catalogue Number	Volume [µl]	Accuracy	Imprecision	Colour Code
613-2891	0.2 - 2	±12.0 - ±1.5 %	≤6.0 - ≤0.7 %	Red
613-7102	0.1 - 2	±40.0 - ±1.5 %	≤12.0 - ≤0.7 %	Red
613-2892	0.5 - 10	±4.0 - ±0.5 %	≤2.8 - ≤0.4 %	Red
613-7103	0.5 - 10	±4.0 - ±0.5 %	≤2.8 - ≤0.4 %	Red
613-2893	2 - 20	±3.0 - ±0.8 %	≤1.5 - ≤0.3 %	Yellow
613-7104	2 - 20	±3.0 - ±0.8 %	≤1.5 - ≤0.3 %	Yellow
613-2894	10 - 100	±1.6 - ±0.8 %	≤0.8 - ≤0.2 %	Yellow
613-7105	10 - 100	±1.6 - ±0.8 %	≤0.8 - ≤0.2 %	Yellow
613-7106	20 - 200	±1.2 - ±0.6 %	≤0.6 - ≤0.2 %	Yellow
613-2896	100 - 1000	±1.6 - ±0.6 %	≤0.4 - ≤0.15 %	Blue
613-7107	100 - 1000	±1.6 - ±0.6 %	≤0.4 - ≤0.15 %	Blue
613-2897	500 - 5000	±1.2 - ±0.5 %	≤0.5 - ≤0.15 %	White
613-7108	500 - 5000	±1.2 - ±0.5 %	≤0.5 - ≤0.15 %	White
613-2898	1000 - 10000	±2.5 - ±0.5 %	≤0.6 - ≤0.2 %	White
613-7109	1000 - 10000	±2.5 - ±0.5 %	≤0.6 - ≤0.2 %	White

Starter Kits

Catalogue Number	Description
613-7110	3 Pack incl. 3 pipettes (2 - 20 / 20 - 200 / 100 - 1000 µl), 2 racks of tips (200 and 1000 µl)
613-6969	VWR Standard 3 Pack (2 - 20 / 20 - 200 / 100 - 1000 µl), without pipette tip samples

Quality System Compliance

Products are manufactured under the **ISO 9001 and ISO 8655** standard. Products are inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable QC requirements.

Product Specifications

Material	Handle part: Polypropylene (PP), Polycarbonate (PC), Polyacetal (POM), Fluoropolymer elastomer, Stainless steel, Brass. Shaft and piston area parts: Polypropylene (PP), Polyvinylidene fluoride (PVDF), Polyphenylene sulfide (PPS), Polyamide (PA), Fluoropolymer elastomer, Stainless steel.
Sterilization	Not Applicable.
ATP Assay	Not Applicable.
DNase & RNase Free	Not Applicable.
BSE/TSE	Not Applicable.
Cytotoxicity	Not Applicable.
Non-Pyrogenic Statement	Not Applicable.
Latex Statement	Not Applicable.
BPA Statement	Not Applicable.

DEHP Statement	Not Applicable.
RoHS	Not Applicable.
REACH Statement	No substances as listed on >Substance of Very High Concern< candidate list published 2023-06-14 are used in manufacturing of the raw materials used to produce this product. No routinely analyse is performed.
Storage Conditions	-25°C to 55°C.
Shelf Life	Not Applicable.

Disclaimer: VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. The purpose of the product is for use in laboratory only.
This document has been produced electronically and is valid without a signature.

Version History

Sl.no	Version number	Description	Month & Year
1	1.0	Created General CoQ using new layout	09-2023