

Material	24549.554
Material description	n-Heptane
Grade	GPR RECTAPUR
CAS Number	142-82-5
Molecular formula	$\text{H}_3\text{C}(\text{CH}_2)_5\text{CH}_3$
Molecular mass	100.20

Characteristics	Specifications
Assay	$\geq 99.0 \%$
IR Spectrum	Passes test
Acidity	$\leq 0.0008 \text{ meq/g}$
Boiling point	$98 - 99^\circ\text{C}$
Density (20/4)	$0.681 - 0.685$
Substances coloured by H_2SO_4	$\leq 90 \text{ APHA}$
Evaporation residue	$\leq 20 \text{ ppm}$

Signature

We certify that this batch conforms to the specifications listed above.

This document has been produced electronically and is valid without a signature.

Isabelle Habay, Head of Laboratory - Briare
VWR International S.A.S.; Z.I. de Vaugereau; FR-45250
Briare; France

For Professional use in Laboratory or Manufacturing. Not for use as an Active Pharmaceutical Ingredient or Food or Animal Feed. Suitability and intended use of the product remains the responsibility of the user.