



M210-LAB Handheld Label Printer

www.vwr.com

Where
Pens Fail.
Printers
Succeed.



When it comes to in-the-lab labeling, illegible or smeared writing is not an option. With the powerful and affordable M210-LAB Handheld Label Printer on hand, you can rely on crisp, long-lasting identification. Once you experience the built-in symbols, durable labels and vial/tube auto-sizing capabilities of the M210-LAB printer, you'll never use a pen again.

Labels Tested for Labs

- Labels and inks built to withstand harsh lab environments, including autoclaves, freezers, liquid nitrogen, and isopropyl alcohol and methyl ethyl ketone chemicals
- Unlike handwritten identification, printed labels reduce errors with easy-to-read text that fits more information in a small space
- Labels engineered for the lab means you never lose your data due to a failed label
- Technical Data Sheets (TDS) available with the specific properties and test results of each material, available at BradyID.com/techdata

Convenient Printer Features

- Vial label shortcut key automatically cuts labels to fit to vial and tube sizes
- Greek symbols and time and date stamps for enhanced sample tracking
- Built-in serialization and bar code symbologies provide more information in a smaller footprint
- Drop-lock-and-print cartridges for quick printing with no setup required
- Easy to use and easy to train others to use
- Variety of power options (use AA batteries, rechargeable battery or AC power supply — all sold separately)

Easy to Own and Maintain

- Label supplies and printer available through top lab distributors
- Low price point allows each technician or lab station to have a printer
- Two-year warranty to ensure it remains in top working condition



M210-LAB Printer

VWR Catalog #	Description
76537-110	M210-LAB Handheld Label Printer

For more information call 1-800-932-5000 or visit www.VWR.com.

Label Materials that Last

Ultra-thin Polyester (B-461)



- Ultra-thin label easily slides in and out of racks and centrifuges
- Ideal for tubes, vials and conicals
- Lasts in liquid nitrogen, autoclaves and hot water baths

VWR Catalog #	Color	Width Inch (mm)	Length Ft (m)
89209-524	Black on White	0.375 (9.53)	21 ft (6.40 m)
89209-528	Black on White	0.750 (19.05)	21 ft (6.40 m)

Clear Polyester Material (B-430)



- For flat or large curved surfaces
- Ideal for general purpose bottles, flasks and dishes
- Lasts in freezers and autoclaves

VWR Catalog #	Color	Width Inch (mm)	Length Ft (m)
89204-988	Black on Clear	0.375 (9.53)	21 ft (6.40 m)
89204-994	Black on Clear	0.500 (12.70)	21 ft (6.40 m)
89205-000	Black on Clear	0.750 (19.05)	21 ft (6.40 m)

* Cartridge is only compatible with the BMP@21-LAB printer.

Polypropylene Label (B-7425)



- Versatile, durable and economical
- Ideal for vials, tubes and other specimen containers
- Lasts in liquid nitrogen

VWR Catalog #	Color	Width Inch (mm)	Length Ft (m)
89233-533	Black on White	0.375 (9.53)	21 ft (6.40 m)
89233-535	Black on White	0.500 (12.70)	21 ft (6.40 m)
89233-537	Black on White	0.750 (19.05)	21 ft (6.40 m)

Durable Vinyl Material (B-595)



- Ideal for equipment, shelves, cabinets and general facility ID
- Available in a variety of colors
- For uneven and rough surfaces

VWR Catalog #	Color	Width Inch (mm)	Length Ft (m)
10067-712	Black on White	0.250 (6.35)	21 ft (6.40 m)
89205-022	Black on White	0.750 (19.05)	21 ft (6.40 m)
10067-714	Black on Yellow	0.250 (6.35)	21 ft (6.40 m)
490006-179	Black on Yellow	0.375 (9.53)	21 ft (6.40 m)
89217-606	Black on Yellow	0.500 (12.70)	21 ft (6.40 m)
89217-608	Black on Yellow	0.750 (19.05)	21 ft (6.40 m)
490006-175	Black on Orange	0.375 (9.53)	21 ft (6.40 m)
10067-706	White on Black	0.375 (9.53)	21 ft (6.40 m)
10067-708	White on Black	0.500 (12.70)	21 ft (6.40 m)
10067-710	White on Black	0.750 (19.05)	21 ft (6.40 m)
490006-177	White on Red	0.375 (9.53)	21 ft (6.40 m)
89205-018	White on Red	0.500 (12.70)	21 ft (6.40 m)
89205-028	White on Red	0.750 (19.05)	21 ft (6.40 m)
490006-171	White on Blue	0.375 (9.53)	21 ft (6.40 m)
89205-014	White on Blue	0.500 (12.70)	21 ft (6.40 m)
89205-024	White on Blue	0.750 (19.05)	21 ft (6.40 m)
10067-634	White on Brown	0.375 (9.53)	21 ft (6.40 m)
10067-636	White on Brown	0.500 (12.70)	21 ft (6.40 m)
10067-638	White on Brown	0.750 (19.05)	21 ft (6.40 m)
490006-173	White on Green	0.375 (9.53)	21 ft (6.40 m)
89205-016	White on Green	0.500 (12.70)	21 ft (6.40 m)
89205-026	White on Green	0.750 (19.05)	21 ft (6.40 m)
10067-629	White on Purple	0.375 (9.53)	21 ft (6.40 m)
10067-630	White on Purple	0.500 (12.70)	21 ft (6.40 m)
10067-632	White on Purple	0.750 (19.05)	21 ft (6.40 m)
10067-700	Black on Grey	0.375 (9.53)	21 ft (6.40 m)
10067-702	Black on Grey	0.500 (12.70)	21 ft (6.40 m)
10067-704	Black on Grey	0.750 (19.05)	21 ft (6.40 m)

* Cartridge is only compatible with the BMP@21-LAB printer.

Chemical-resistant Polyester Label (B-488)



- Ideal for slides
- Lasts in freezers and autoclaves, and withstands DMSO, ethanol and short-term xylene exposure

VWR Catalog #	Color	Width Inch (mm)	Length Ft (m)
89209-530	Black on White	0.375 (9.53)	21 ft (6.40 m)
89209-532	Black on White	0.500 (12.70)	21 ft (6.40 m)
89209-534	Black on White	0.750 (19.05)	21 ft (6.40 m)

Nylon Cloth Label (B-499)



- Easy-to-handle label for non-rigid and curved surfaces
- Ideal for bags, bottles and tubes
- Lasts in liquid nitrogen, freezers and autoclaves

VWR Catalog #	Color	Width Inch (mm)	Length Ft (m)
89204-986	Black on White	0.375 (9.53)	16 ft (4.88 m)
89204-992	Black on White	0.500 (12.70)	16 ft (4.88 m)
89204-998	Black on White	0.750 (19.05)	16 ft (4.88 m)