


Physics

From classic and cutting-edge apparatus and demos, to time-saving activities, durable equipment, and essential supplies for your lab - we have what you need to turn physics lessons into physics connections.



Specifications:

DC Power Supply 20V,5A

- DC output: 0-20 V, 0-5 A
- Output power: 100 W
- Stability under full load $\leq 0,01\% + 5 \text{ mV}$, $\leq 0,2\% + 5 \text{ mA}$
- Residual ripple $\leq 1 \text{ mV}$, 3 mA
- Display: 2 x 3 digit LED
- Terminals: 4 mm safety sockets
- Dimensions: approx. 130x150x300 mm³
- Weight: approx. 4.7 kg