

Certificate of Conformity

VWR Catalogue Number	See List
Description	VWR® Forceps
Country of Origin	Manufactured in United Kingdom
Date of Issue (yyyy-mm-dd)	2026-05-25

Catalogue Number	Type	Length [mm]
232-0150	Square end	100
232-0151	Sharp end	100
232-0152	Square end	150
232-0153	Sharp end	150
232-0154	Square end	200
232-0155	Sharp end	200

Quality System Compliance

This certificate confirms that the product has been sourced from a supplier or manufacturer who has declared their Quality Management System (QMS) complies with **ISO 9001:2015**.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable QC requirements.

Product Specifications

Material	Virgin Polytetrafluoroethylene (PTFE). Raw Material complies with Directive FDA 21 CFR 177.1550 and complies with USP Class VI plastic (121 °C).
Temperature Resistance	–200 °C to +280 °C.
Sterilization	Non-Sterile.
ATP Assay	Not Applicable

DNase & RNase Free	Not Applicable
BSE/TSE	No use of any raw material produced from or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE).
Cytotoxicity	Not Applicable
Non-Pyrogenic Statement	Not Applicable
Latex Statement	This product is not made from natural rubber latex.
BPA Statement	Bisphenols are not used in the manufacture of the raw material and are not expected to be present.
DEHP Statement	The substance Phthalates are not used or intentionally incorporated in the manufacture of the raw materials or final product.
RoHS	Product does not contain any of the restricted substance above the maximum concentration values as listed in Annex II of RoHS 2 Directive 2011/65/EU, as amended by EU 2015/863.
REACH Statement	The material used for the production of this container does not contain more than 0.1 % (w/w) of any substance of very high concern listed in the ECHA candidate list.
Storage Conditions	No special storage conditions required.
Shelf Life	Indefinite when stored under ideal conditions.

Disclaimer: VWR's Quality Management System (QMS) is ISO 9001:2015 certified as a wholesale distributor of finished products, ensuring effective oversight of finished product handling, traceability, and distribution processes. VWR is not a manufacturer and does not engage in manufacturing, co-manufacturing or processing activities such as design, packaging, labeling, testing or quality release. Consequently, the supplier or manufacturer's QMS certification may differ from VWR's.

VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. This product is intended for research use only. It is not for use in diagnostic procedures, therapeutic applications, or any human or clinical use.

This document has been produced electronically and is valid without a signature.

Version History

Sl.no	Version number	Description	Month & Year
1	1.0	Created General CoC using new layout	05-2026