


Physics

From classic and cutting-edge apparatus and demos, to time-saving activities, durable equipment, and essential supplies for your lab - we have what you need to turn physics lessons into physics connections.



Specifications:

DC Power Supply 0-500V

- 500 V output:
 - Voltage: 0 – 500 V DC, max. 50 mA
 - Stability at full load: $\leq 0.01\% \pm 100 \text{ mV}$
 - Residual ripple: $\leq 20 \text{ mV}$
- 50 V output:
 - Voltage: 0 – 50 V DC, max. 50 mA
 - Stability at full load: $\leq 0.1\% \pm 30 \text{ mV}$
 - Residual ripple: $\leq 50 \text{ mV}$
- 12 V output:
 - Voltage: 0 – -12 V DC, max. 4 A
 - Stability at full load: $\leq 0.1\% \pm 30 \text{ mV}$
- 8 V output:
 - Voltage: 0 – 8 V DC, max. 3 A
 - Stability at full load: $\leq 0.1\% \pm 30 \text{ mV}$
- Displays: Analogue, class 2
- Connections: 4 mm safety sockets
- Power consumption: 50 VA