



Specification

1.00898.0010 E.coli/Coliform Selective-Supplement

Specification

Supplement tested in Chromocult Coliform agar art. no. 1.10426.

Recovery on test medium (Escherichia coli DSM 502)	≥ 70	%
--	------	---

Recovery on test medium (Citrobacter freundii ATCC 8090)	≥ 70	%
--	------	---

Growth (Aeromonas hydrophila ATCC 7966 (WDCM 00063))	inhibited
--	-----------

Growth (Pseudomonas aeruginosa ATCC 27853 (WDCM 00025))	inhibited
---	-----------

Growth (Bacillus licheniformis ATCC 14580 (WDCM 00068))	inhibited
---	-----------

Growth (Staphylococcus aureus ATCC 6538 (WDCM 00032))	inhibited
---	-----------

Colony color (Escherichia coli DSM 502)	blue
---	------

Colony color (Citrobacter freundii ATCC 8090)	red
---	-----

Incubation: 24 hrs.; 35 °C; aerobic.

Stefanie Fischer

Responsible laboratory manager quality control

This document has been produced electronically and is valid without a signature.