



Product Specification Sheet

VWR MICROSCOPE SLIDES, CUT EDGES, FLOAT GLASS

VWR International certifies that the specified products meet all VWR specifications and quality criteria for release. The products comply in every respect with the relevant particular specifications, drawings, relevant standards and regulations in force.

VWR International certifie que les produits mentionnés respectent toutes les spécifications et critères de qualité VWR. Les produits satisfont aux spécifications particulières, aux plans, aux normes et règlements en vigueur s'y rapportant.

VWR International bestätigt hiermit, daß die angegebenen Produkte alle VWR Spezifikationen und Qualitätskriterien zur Freigabe erfüllen. Die Produkte entsprechen in jeder Hinsicht den maßgeblichen Spezifikationen, den Zeichnungen sowie den relevanten Normen und gültigen Vorschriften.

Cat. No.	Description	Pk
631-1550	Plain	50
631-1551	Frosted	50

Material: Transparent Float glass

WxDxH: 76x26x1 mm

Other characteristics:

- Chemical durability, resistance to the atmosphere
- Flatness should be no more than 0.05 mm
- The thickness of the slides should be symmetrical and tolerance of thickness variation should be more than 0.05 mm
- The width of the frosted end should be 19mm±2.0 mm

Disclaimer: VWR wants to state that these specifications can only be regarded as general information and will not discharge the user from his obligation to ensure the product is suitable for the intended use - If there is any doubt ask us.

Date of issue: **April 2013**