

Material	83612.360
Material description	Water
Grade	GPR RECTAPUR produced with reversed osmosis and CDI
CAS Number	7732-18-5
Molecular formula	H ₂ O
Molecular mass	18.02

Characteristics	Specifications
Acidity or alkalinity	Passes test
Oxidisable substances	Passes test
Conductivity (20°C)	≤ 1.5 µS/cm
Evaporation residue	≤ 1 ppm
Heavy metals (as Pb)	≤ 0.1 ppm
Cl (Chloride)	≤ 0.1 ppm
NH ₄ (Ammonium)	≤ 0.2 ppm
NO ₃ (Nitrate)	≤ 0.2 ppm
SO ₄ (Sulphate)	≤ 1 ppm
Ca (Calcium)	≤ 0.1 ppm
Na (Sodium)	≤ 0.1 ppm

Signature

We certify that this batch conforms to the specifications listed above.

This document has been produced electronically and is valid without a signature.

Jessie Pijck, Head of Laboratory - Haasrode
VWR International bv; Geldenaaksebaan 464; BE-3001
Leuven; Belgium

For Professional use in Laboratory or Manufacturing. Not for use as an Active Pharmaceutical Ingredient or Food or Animal Feed. Suitability and intended use of the product remains the responsibility of the user.