

# Certificate of Conformity

|                                   |                            |
|-----------------------------------|----------------------------|
| <b>VWR Catalogue Number</b>       | See List                   |
| <b>Description</b>                | Pipette Tip, Reload System |
| <b>Country of Origin</b>          | Manufactured in USA        |
| <b>Date of Issue (yyyy-mm-dd)</b> | 2021-08-17                 |

| Catalogue Number | Tip Style                    | Sterility   | Volume<br>[ µl ] |
|------------------|------------------------------|-------------|------------------|
| 613-1068         | Micro, graduated             | Non-Sterile | 0.1 - 10         |
| 613-1069         | Micro, extra-long, graduated | Non-Sterile | 0.1 - 10         |
| 613-1067         | Micro, Eppendorf type        | Non-Sterile | 0.5 - 20         |
| 613-1065         | Graduated, bevelled          | Non-Sterile | 1 - 200          |
| 613-1064         | Graduated, bevelled (yellow) | Non-Sterile | 1 - 200          |
| 613-1066         | Micro                        | Non-Sterile | 1 - 200          |
| 613-1063         | Graduated, bevelled          | Non-Sterile | 5 - 300          |
| 613-1062         | Graduated                    | Non-Sterile | 100 - 1250       |
| 613-1075         | Micro, graduated             | Non-Sterile | 0.1 - 10         |
| 613-1076         | Micro, extra-long, graduated | Non-Sterile | 0.1 - 10         |
| 613-1074         | Micro, Eppendorf type        | Non-Sterile | 0.1 - 20         |
| 613-1072         | Graduated, bevelled          | Non-Sterile | 1 - 200          |
| 613-1073         | Micro                        | Non-Sterile | 1 - 200          |
| 613-1071         | Graduated, bevelled          | Non-Sterile | 5 - 300          |
| 613-1070         | Graduated                    | Non-Sterile | 100 - 1250       |
| 613-1082         | Micro, graduated             | Sterile     | 0.1 - 10         |
| 613-1083         | Micro, extra-long, graduated | Sterile     | 0.1 - 10         |
| 613-1081         | Micro, Eppendorf type        | Sterile     | 0.5 - 20         |
| 613-1079         | Graduated, bevelled          | Sterile     | 1 - 200          |
| 613-1080         | Micro                        | Sterile     | 1 - 200          |
| 613-1078         | Graduated, bevelled          | Sterile     | 5 - 300          |
| 613-1077         | Graduated                    | Sterile     | 100 - 1250       |

| Catalogue Number | Tip Style           | Sterility   | Volume<br>[ µl ] |
|------------------|---------------------|-------------|------------------|
| 613-1093         | Micro, graduated    | Non-Sterile | 0.1 - 10         |
| 613-1090         | Graduated           | Non-Sterile | 1 - 200          |
| 613-1089         | Graduated (yellow)  | Non-Sterile | 1 - 200          |
| 613-1091         | Micro               | Non-Sterile | 1 - 200          |
| 613-1092         | Micro (yellow)      | Non-Sterile | 1 - 200          |
| 613-1088         | Graduated, bevelled | Non-Sterile | 5 - 300          |
| 613-1087         | Graduated           | Non-Sterile | 100 - 1250       |
| 613-1100         | Micro, graduated    | Sterile     | 0.1 - 10         |
| 613-1097         | Graduated           | Sterile     | 1 - 200          |
| 613-1096         | Graduated (yellow)  | Sterile     | 1 - 200          |
| 613-1098         | Micro               | Sterile     | 1 - 200          |
| 613-1099         | Micro (yellow)      | Sterile     | 1 - 200          |
| 613-1095         | Graduated, bevelled | Sterile     | 5 - 300          |
| 613-1086         | Empty rack          | -           | 10 - 20          |
| 613-1085         | Empty rack          | -           | 200 - 300        |
| 613-1084         | Empty rack          | -           | 1000 - 1250      |

### Quality System Compliance

Products are manufactured under the **ISO 9001 & ISO 13485** standard. Products are inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP.

### QC Testing

Representative products samples are collected and inspected in accordance with current applicable QC requirements.

**PCR Inhibitor:** Products tested to be considered free of PCR inhibitors.

## Product Specifications

|                                |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                     |
|--------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Material</b>                | High Grade Polypropylene.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| <b>Sterilization</b>           | Product labelled as sterile are Electron Beam irradiated according to guidelines of ANSI/AAMI/ISO 11137:2006 and meets SAL 10 <sup>-6</sup> with a specified dose range of min. 15.8 kGy and max. 40.0 kGy which is not applicable for Non-sterile products.                                                                                                                                                                                                                                                                                                                                        |
| <b>ATP Assay</b>               | This product has been tested and is ATP free. The assay detection limit is <10 <sup>-13</sup> mg/μL.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                |
| <b>DNase &amp; RNase Free</b>  | This product is free of any detectable DNase/RNase contamination. The assay detection limit is <10 <sup>-7</sup> Units/μL for DNase and <10 <sup>-9</sup> Units/μL for RNase.                                                                                                                                                                                                                                                                                                                                                                                                                       |
| <b>BSE/TSE</b>                 | <p>No use of any raw material produced from or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE).</p> <p>Tallow derived additives may be used in the manufacturing of the raw material. The material used in the product comply with the Regulations 1069/2009/EC, and 142/2011/EC and the "Note for guidance EMEA 410/01, rev.3".</p> |
| <b>Cytotoxicity</b>            | Not applicable.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                     |
| <b>Non-Pyrogenic Statement</b> | The acceptance level for product is <0.06 EU/mL.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
| <b>Latex Statement</b>         | This product is not made from natural rubber latex.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| <b>BPA Statement</b>           | Not applicable.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                     |
| <b>DEHP Statement</b>          | Please see RoHS information below.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| <b>RoHS</b>                    | No substances (Lead, Cadmium, Mercury, Hexavalent Chromium (Cr6+), Poly Brominated Biphenyls (PBB), Poly Brominated Diphenyl ethers (PBDE), Bis(2-Ethyhexyl) Phthalate (DEHP), Benzyl Butyl Phthalate (BBP), Dibutyl Phthalate (DBP), Diisobutyl Phthalate (DIBP)) are used in manufacturing the raw materials and final product. No routinely analyse is performed.                                                                                                                                                                                                                                |

|                           |                                                                                                                                                                                                                  |
|---------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>REACH Statement</b>    | No substances as listed on >Substance of Very High Concern< candidate list published 2021-07-08 are used in manufacturing of the raw materials used to produce this product. No routinely analysis is performed. |
| <b>Storage Conditions</b> | Store in a room temperature.                                                                                                                                                                                     |
| <b>Shelf Life</b>         | 5 years from the date of Manufacture which is not applicable for Non-sterile products.                                                                                                                           |

Disclaimer: VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. The purpose of the product is for use in laboratory only.  
This document has been produced electronically and is valid without a signature.