

Certificate of Quality

VWR Catalogue Number	See List
Description	VWR® Sterile Cell Culture Flasks
Country of Origin	Manufactured in China
Date of Issue (yyyy-mm-dd)	2021-03-16

VWR Catalogue Number	Capacity [mL]	Cap Type	Culture Area [cm²]
10062-862	250	Standard	75
10062-860	250	Vented	75
10062-870	25	Standard	12.5
10062-868	25	Vented	12.5
10062-874	50	Standard	25
10062-872	50	Vented	25
10062-866	600	Standard	182
10062-864	600	Vented	182
10062-886	850	Standard	300
10062-884	850	Vented	300

Quality System Compliance

Products are manufactured under the **ISO 9001:2015 & ISO 13485:2016** standard. Products are inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable QC requirements.

Product Specifications

This product is Class I category device as defined by the FDA in 21CFR Parts 862-892.

Material	Flask	Polystyrene (PS)
	Cap	High density PolyEthylene (HDPE)

Sterilization	Products are EB (Electron Beam) irradiated and dose released upon ISO 11137 recommended practices in effect at the time of validation and they meet a minimum requirement of 10 ⁻⁶ SAL with a specified dose range of 15 - 30 kGy.
ATP Assay	Not applicable.
DNase & RNase Free	This product is free of any detectable DNase/RNase contamination.
BSE/TSE	No use of any raw material produced from or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE).
Cytotoxicity	Testing is conducted to quality all material resins using ISO 10993 standards for cytotoxicity and have been shown to be non-toxic.
Non-Pyrogenic Statement	The acceptance level for product is 0.5 EU/ml or less than 20 EU/device. (TAL Gel Clot Method).
Latex Statement	The product is latex free.
BPA Statement	Bisphenol are not used in the manufacture of the raw material and are not expected to be present.
DEHP Statement	Not applicable.
RoHS	No substances (Lead, Cadmium, Mercury, Hexavalent Chromium (Cr ⁶⁺), Poly Brominated Biphenyls (PBB), Poly Brominated Diphenyl ethers (PBDE), Benzyl Butyl Phthalate (BBP), Dibutyl Phthalate (DBP), Diisobutyl Phthalate (DIBP), DEHP (di(2- ethylhexyl) phthalate)) are used in manufacturing the raw materials and final product. No routinely analyse is performed.
REACH Statement	Not applicable.
Storage Conditions	Store at room temperature.
Shelf Life	3 Years.

Disclaimer: VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. The purpose of the product is for use in laboratory only.
This document has been produced electronically and is valid without a signature.