

Certificate of Conformity

VWR Catalogue Number	216-0892
Description	Conical Faeces / specimen containers, with screw cap, Translucent, 120 ml Capacity, 71.8 mm Height, 56 mm External Diameter, Red screw cap and spatula, Sterile
Lot Number	21C0990
Date of Manufacture (yyyy-mm)	2021-03
Date of Sterilization (yyyy-mm-dd)	2021-03-10
Sterilization Expiry Date (yyyy-mm-dd)	2026-03-10
Country of Origin	Manufactured in Italy
Date of Issue (yyyy-mm-dd)	2021-12-16

Quality System Compliance

Products are manufactured under the **ISO 9001:2015**, **ISO 13485:2016** standard. Products are inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable QC requirements.

Product Specifications

These products are in compliance with CE IVD Directive 98/79/EC.

Material	Polypropylene.
Sterilization	Products are sterilized using Beta rays and meets SAL 10 ⁻⁶ with a specified dosage range of 17 – 20 kG.
ATP Assay	Not applicable.
DNase & RNase Free	Not applicable.

BSE/TSE	No use of any raw material produced from or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE).
GMO Statement	Genetically modified materials (GMO) have not been intentionally incorporated in the manufacturing of the product.
Cytotoxicity	Not applicable.
Non-Pyrogenic Statement	Not applicable.
Latex Statement	The products are not made from natural rubber latex.
BPA Statement	Bisphenol A (BPA) has not been intentionally incorporated in the manufacturing of the products.
DEHP Statement	DEHP has not been intentionally incorporated in the manufacturing of the products.
RoHS	Not applicable.
REACH Statement	No substances as listed on >Substance of Very High Concern< candidate list published 2021-07-08 are used in manufacturing of the raw materials used to produce this product. No routinely analysis is performed.
Storage Conditions	Store in a dry place, away from sources of heat and direct sunlight.

Disclaimer: VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. The purpose of the product is for use in laboratory only.
This document has been produced electronically and is valid without a signature