



VWR® SPEC-WIPE POLYESTER

As per the manufacturer, the below product meets the following criteria:

VWR Catalog Number:	76377-480
Lot Number:	798376-0-0-1
Description:	VWR SPEC-WIPE POLYESTER 12X12 IN CS300
Date of Manufacture:	August 29, 2020
Expiration Date:	August 2022
Sterilization Lot Number:	58567A
Specified Dosage (kGy):	Max: 50.0; Min: 25.0
Delivered Dosage (kGy):	Max: 42.8; Min: 27.0
Country of Origin:	US

Quality System Compliance

The above product was manufactured in accordance with prescribed procedures, specifications and under ISO 9001 standard.

Product Specifications:

Sterilization Process:	The above product has been dosimetrically released as sterile. All aspects of the manufacturing process are within specifications and the specified minimum dose has been delivered to the product. This is in accordance with procedures outlined in ANSI/AAMI/ISO 11137-2 Method VDMax: <i>Sterilization of healthcare products – Requirements for validation and routine control – Radiation Sterilization</i> , to a Sterility Assurance Level of 10^{-6} .
BSE/TSE Statement:	This product is manufactured wholly from synthetic materials and do not contain any raw materials produced from or substances derived of animal origin. The manufacturing process does not use any ingredient of animal origin, nor do the materials come into contact with animal products during storage and transportation. The products are free from Transmissible Spongiform Encephalopathy (TSE) and Bovine Spongiform Encephalopathy (BSE).
Non-Pyrogenic:	N/A
DNase & RNase Free:	N/A
Cytotoxicity:	N/A
Latex Free:	Latex or latex derivatives are not used in the production of this product.

Signed:



Jamie Ethier
VP Global Quality
VWR, Part of Avantor

Date: October 12, 2020