

# Certificate of Conformity

|                                   |                                                                    |
|-----------------------------------|--------------------------------------------------------------------|
| <b>VWR Catalogue Number</b>       | See List                                                           |
| <b>Description</b>                | VWR® Spec-Wipe® 4, Cleanroom Wipes,<br>229 mm length, 229 mm width |
| <b>Country of Origin</b>          | Manufactured in China                                              |
| <b>Date of Issue (yyyy-mm-dd)</b> | 2026-04-15                                                         |

| Catalogue Number | Description                                               | Packed      |
|------------------|-----------------------------------------------------------|-------------|
| 115-0036         | Spec-Wipe® 4, dry                                         | 150 per bag |
| 115-0039         | Spec-Wipe® 4, saturated with 70% IPA, 30% deionised water | 30 per bag  |

## Quality System Compliance

This certificate confirms that the product has been sourced from a supplier or manufacturer who has declared their Quality Management System (QMS) complies with **ISO 9001:2015**.

## QC Testing

Representative products samples are collected and inspected in accordance with current applicable QC requirements.

## Product Specifications

|                               |                 |
|-------------------------------|-----------------|
| <b>Material</b>               | 100% polyester. |
| <b>Sterilization</b>          | Non-Sterile.    |
| <b>ATP Assay</b>              | Not Applicable. |
| <b>DNase &amp; RNase Free</b> | Not Applicable. |

|                                |                                                                                                                                                                                                                                                                                                                                                             |
|--------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>BSE/TSE</b>                 | No use of any raw material produced from or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE). |
| <b>Cytotoxicity</b>            | Not Applicable.                                                                                                                                                                                                                                                                                                                                             |
| <b>Non-Pyrogenic Statement</b> | Not Applicable.                                                                                                                                                                                                                                                                                                                                             |
| <b>Latex Statement</b>         | Latex or latex derivatives are not used in the production of this product.                                                                                                                                                                                                                                                                                  |
| <b>BPA Statement</b>           | Not Applicable.                                                                                                                                                                                                                                                                                                                                             |
| <b>DEHP Statement</b>          | Not Applicable.                                                                                                                                                                                                                                                                                                                                             |
| <b>RoHS</b>                    | Not Applicable.                                                                                                                                                                                                                                                                                                                                             |
| <b>REACH Statement</b>         | Not Applicable.                                                                                                                                                                                                                                                                                                                                             |
| <b>Storage Conditions</b>      | Stored in standard warehouse conditions within an environment of ambient temperature and humidity.                                                                                                                                                                                                                                                          |
| <b>Shelf Life</b>              | No shelf life.                                                                                                                                                                                                                                                                                                                                              |

Disclaimer: VWR's Quality Management System (QMS) is ISO 9001:2015 certified as a wholesale distributor of finished products, ensuring effective oversight of finished product handling, traceability, and distribution processes. VWR is not a manufacturer and does not engage in manufacturing, co-manufacturing or processing activities such as design, packaging, labeling, testing or quality release. Consequently, the supplier or manufacturer's QMS certification may differ from VWR's.

VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. This product is intended for research use only. It is not for use in diagnostic procedures, therapeutic applications, or any human or clinical use.

This document has been produced electronically and is valid without a signature.

## Annex

| Attribute; (units)                                          | Value ** | Test Method                 |
|-------------------------------------------------------------|----------|-----------------------------|
| Basis weight; (g/m <sup>2</sup> )                           | 136      |                             |
| Absorbency in water                                         |          |                             |
| Intrinsic; (mL/g)                                           | 2.3      | IENT-RP-CC004.2, Sec. 7.1   |
| Extrinsic; (mL/m <sup>2</sup> )                             | 312      | IENT-RP-CC004.2, Sec. 7.1   |
| Sorptive rate; (seconds)                                    | <1.5     |                             |
| Extractables                                                |          | IENT-RP-CC004.2, Sec. 6.1.2 |
| In deionized water; (g/m <sup>2</sup> )                     | 0.042    |                             |
| In isopropyl alcohol; (g/m <sup>2</sup> )                   | 0.015    |                             |
| Specific ions                                               |          | IENT-RP-CC004.2, Sec. 6.2.2 |
| Sodium; (ppm)                                               | 0.39     |                             |
| Chloride (ppm)                                              | 0.45     |                             |
| Particle generation                                         |          | IENT-RP-CC004.2, Sec. 5.1   |
| P <sub>0</sub> ≥ 0.5µm; (x10 <sup>6</sup> /m <sup>2</sup> ) | 4.8      |                             |
| Fibers > 100µm; (x10 <sup>3</sup> /m <sup>2</sup> )         | 3.3      | IENT-RP-CC004.2, Sec. 5.2   |

\*\* ND = None detected; levels are below detection limit of test equipment

## Version History

| Sl.no | Version number | Description                          | Month & Year |
|-------|----------------|--------------------------------------|--------------|
| 1     | 1.0            | Created General CoQ using new layout | 04-2026      |
|       |                |                                      |              |
|       |                |                                      |              |