



Product Specification Sheet

VWR TAPE, POLYETHYLENE, INDUSTRIAL GRADE

VWR International certifies that the specified products meet all VWR specifications and quality criteria for release. The products comply in every respect with the relevant particular specifications, drawings, relevant standards and regulations in force.

VWR International certifie que les produits mentionnés respectent toutes les spécifications et critères de qualité VWR. Les produits satisfont aux spécifications particulières, aux plans, aux normes et règlements en vigueur s'y rapportant.

VWR International bestätigt hiermit, daß die angegebenen Produkte alle VWR Spezifikationen und Qualitätskriterien zur Freigabe erfüllen. Die Produkte entsprechen in jeder Hinsicht den maßgeblichen Spezifikationen, den Zeichnungen sowie den relevanten Normen und gültigen Vorschriften.

Cat. No.	Width	Length
817-0100	25.5 mm	32.9 m
817-0095	25.5 mm	32.9 m
817-0096	25.5 mm	32.9 m
817-0099	25.5 mm	32.9 m
817-0097	25.5 mm	32.9 m
817-0098	25.5 mm	32.9 m
817-0101	25.5 mm	32.9 m

Other Specifications:

Backing: **White, pigmented and transparent low-density polyethylene film**

Adhesive: **Acrylic based pressure sensitive**

Adhesion to steel: **56 oz./inch width**

Tensile Strength: **12 lbs./inch width**

Total thickness: **0.0055"**

Elongation at break: **400 %**

Optimum Temp range: **-40°F to 190°F/ -40°C to 88°C**

VWR International bvba/sprl, Researchpark Haasrode 2020, Geldenaaksebaan 464, 3001 Leuven, Belgium

<http://www.vwr.com>





Product Specification Sheet

VWR TAPE, POLYETHYLENE, INDUSTRIAL GRADE

Abrasion Resistance: **Good**

Available colors: **White, Black, Orange, Yellow, Red, Transparent, Blue and Green**

General Information:

A low-density polyethylene tape low in halogens and leachable chlorides. Processed and packaged in a controlled environment under laminar flow. Ideal for packaging and transportation of wafer boxes as the polyethylene film and acrylic adhesive minimize outgassing and offer good clean peel properties. This tape has excellent chemical resistance and good conformability.

Disclaimer: VWR wants to state that these specifications can only be regarded as general information and will not discharge the user from his obligation to ensure the product is suitable for the intended use
- If there is any doubt ask us.

Date of issue: **October 2012**

VWR International bvba/sprl, Researchpark Haasrode 2020, Geldenaaksebaan 464, 3001 Leuven, Belgium

<http://www.vwr.com>



Version: 2.0